User Manual

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Manuale d'uso

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VIXWIN™
PLATINUM
with VxDICOM Module

Imaging Software

GENDEX®
Imaging Excellence Since 1893

Warranty Statement

There are no warranties, express or implied, with respect to the contents of this document, and all information provided herein is provided "as is". Gendex reserves the right to periodically change the information that is contained in this document; however, Gendex makes no commitment to provide any such changes in a timely manner or at all.

Limitation of Liability

In no event shall Gendex or its subsidiaries or affiliates or parent be liable to any party for any direct, indirect, special, or consequential damages in connection with your use of this document, including, without limitation, loss of business revenue or earnings, lost data, damages caused by delays, lost profits, or a failure to realize expected savings, even if Gendex was expressly advised of the possibility of such damages.

Trademarks

Gendex and related marks, images, and symbols are the exclusive properties of, and trademarks of, Gendex Dental Systems. All other brands and marks are the properties of their respective owners. Imaging Technology provided under License by AccuSoft Corporation. ImageGear 1996 - 2001 by AccuSoft Corporation. All Rights Reserved.

Copyright

Documentation for VixWin Platinum and the operating software are copyrighted with all rights reserved. Under the copyright laws, this documentation may not be reproduced, transmitted, transcribed, or translated into any human or computer language in whole or part without the prior written permission of the copyright holder.

About Gendex

To learn more about Gendex, access the Gendex Dental System Website www.gendex.com

Language

The original language of this manual is English.

Table of Contents

Chapter 1 - Product Description

Intended Use	1-2
--------------------	-----

Chapter 2 - Installation and Configuration

Software Installation	2-1
Install	2-1
Uninstall	2-1
Registration	2-1
Configuration	2-1
Storage SCU	2-2
Print	2-3
Worklist	2-3
Query/Retrieve	2-3
Storage SCP	2-4
Optional Parameter Settings	2-4

Chapter 3 - Operation with Remote Servers

Modality Worklist Server	3-1
Image Storage	3-2
Query/Retrieve	3-3
Patient Information	3-5
DICOM Removable Media	3-5
DICOM Print	3-6

Chapter 4 - Glossary

Applicable Standards	4-2
----------------------------	-----

CHAPTER 1 **Product Description**

VxDICOM is an add-on module for VixWin Platinum.

Two versions of VxDICOM are available: Standard and Advanced.

With either of the two versions you can:

- Store images on a remote DICOM storage server (PACS).
- Print images to a DICOM Printer.

In addition, with VxDICOM Advanced you can:

- Select patients/studies from a remote DICOM Worklist Server (RIS).
- Retrieve images from a remote DICOM storage server (PACS).

VxDICOM operates with these operating systems:

- Microsoft® Windows XP®
- Microsoft® Windows Vista® (Business, 32 bit)
- Microsoft® Windows 7® Professional (32 and 64 bit)
- Microsoft® Windows 8® Professional (32 and 64 bit)

This manual can be viewed by choosing **Help > VxDICOM Manual**.

A DICOM Conformance Statement document for the VxDICOM product can be found in the VixWin Platinum CD-ROM under Manuals\VxDICOM folder.

The CE symbol ensures that the product here in specified meets the provisions of European Council Directive 93/42 EEC concerning medical devices.

VxDICOM uses the OFFIS DICOM Toolkit DCMTK © 1993-2006, OFFIS e.V.

Intended Use

VixWin Platinum is an application program used for dental imaging. VxDICOM adds DICOM communication features to the VixWin Platinum product:

- Modality Worklist.
- Store images to DICOM Server.
- Query/Retrieve images from DICOM Server.
- Print images to DICOM Printer.

CHAPTER 2 **Installation and Configuration**

Software Installation

Install

1. Ensure VixWin Platinum is already installed.
2. Insert VixWin installation CD into CD drive.
3. Setup automatically starts.
4. Click Modules, then VxDICOM.

Uninstall

1. Open the Windows Control Panel.
2. Open Add/Remove Programs.
3. Select VxDICOM from the list.
4. Click the Add/Remove button. Follow the on-screen instructions.

Registration

See the VixWin Platinum User Manual for an explanation of how to purchase the license and register the software. When purchasing the license for VxDICOM, you have a choice to purchase the VxDICOM Advanced or VxDICOM Standard versions.

Configuration

Configure VxDICOM on each workstation. Double click the VxDICOM Config icon on the desktop to display the configuration panel below.

Menu Name	IP Address (*)	Port (**)	Calling Title	Called Title
Storage 1	172.10.20.39	204	VIXWIN	ArchiveSCP
Print 1	localhost	204	VIXWIN	DICOM_Print
Worklist	localhost	106	VIXWIN	DICOM_Worklist
Query/Retrieve 1	172.10.20.39	204	VIXWIN	ArchiveSCP
Storage SCP 1	localhost	250	VIXWIN	VIXWIN

(*) IPv4 Address or Host Name (**) Insert decimal values only

Patient Database: VIXWIN

Image Database: VIXWIN

Windows Layout: VIXWIN

Full Mouth: VIXWIN

Advanced Settings:

- ☐ Use Patient Health Code as DICOM Patient ID
- ☒ Patient Information ☐ Auto Patient Information
- ☐ Enable Log Console ☐ Auto Worklist

Storage Modality IOD	Save Treated Image	Save Annotations
Intra Oral: INTRA-ORAL X-RAY(IQ)	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Panoramic: PANORAMIC(PX)	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Cephalometric: DIGITAL X-RAY(DX)	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Camera: SECONDARY CAPTURE(SC)	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Other: INTRA-ORAL X-RAY(IQ)	<input checked="" type="checkbox"/>	<input type="checkbox"/>

DEFAULT SAVE EXIT

Configure remote DICOM Servers to select patients and store, retrieve and print images. For each DICOM functionality (Storage, Print, Worklist, Query\Retrieve and Storage SCP) enter:

- IP Address (or hostname if DNS is active);
- Port Number;
- Calling Title (Application Entity name);
- Called Title (Application Entity name).

DICOM functionalities are described below.

Storage SCU

Parameters required to interface with a DICOM Server providing PACS Services (up to three can be specified).

- Menu Name: Enter text for the VixWin menu name.
- IP Address: IP address or host name of DICOM Storage Server location.
- Port: Port number.

- Called Title: Storage Application Entity Name (i.e. DICOM_Storage).
- Calling Title: VixWin Application Entity Name (i.e. VixWin).

Print

Parameters required to interface with a DICOM printer (up to three can be specified).

- Menu Name: Enter text for the VixWin menu name.
- IP Address: IP address or host name of the DICOM Printer.
- Port: Port number.
- Called Title: Print Application Entity Name.
- Calling Title: VixWin Application Entity Name (i.e. VixWin).

Worklist

Parameters required to interface with a DICOM Server providing worklist (RIS) services.

- IP Address: IP address or host name of DICOM Worklist Server location.
- Port: Port number.
- Called Title: Worklist Application Entity Name (i.e. "DICOM_Worklist").
- Calling Title: VixWin Application Entity Name (i.e. VixWin).

Query/Retrieve

Parameters required to interface with a DICOM Server providing PACS services (up to three can be specified).

- Menu Name: Enter text for the VixWin menu name.
- IP Address: IP address or host name of DICOM Q/R server location.
- Port: Port number.
- Called Title: Q/R Application Entity Name (i.e. "DICOM_Q/R").
- Calling Title: VixWin Application Entity Name (i.e. VixWin).

Storage SCP

Parameters required to interface with a DICOM Server providing PACS services.

- IP Address: IP address or host name of VixWin location.
- Port: Port number.
- Called Title: Application Entity Name (i.e. VixWin).

NOTE: To find a workstation IP address, open a command window and run “IPCONFIG”.

Optional Parameter Settings

It is also possible to configure the following parameters (defaults are underlined):

Patient Database

VIXWIN: Handle patient selection in the same way as VixWin Platinum (Bridge mode or stand-alone database mode – see the VixWin Platinum User Manual).

DICOM Worklist: Use a DICOM Worklist Server.

Image Database

VIXWIN: Images are saved locally in the VixWin database. Images can still be sent individually to a DICOM Storage Server.

DICOM: Use a DICOM Storage Server.

DICOM+VIXWIN: Images are saved both locally in the VixWin database and also on a DICOM Storage Server.

NOTE: When the local VixWin database is selected, images are saved as:

C:\VxImages\<PATIENTID>nn.ext

Where nn is a progressive number and ext is the file format extension.

Use Patient Health Code as DICOM Patient ID

ON: Use Patient Health Code as Patient ID (stand-alone database mode).

OFF: Use VixWin generated Patient ID (stand-alone database mode).

Patient Information

ON: Allows manual entry of patient information.

OFF: Disable this function.

Enable Log Console

OFF: Hides the DICOM SCU Console during storage operation.

ON: Shows the DICOM SCU Console during storage operation.

Auto Patient Information

ON: Patient Information dialog is displayed each time VixWin is run.

OFF: Not displayed.

AutoWorklist

ON: The Worklist window is opened each time VixWin is run.

OFF: Not opened.

Storage Modality IOD

It is possible to specify the image saving modality for different image types:

- intraoral X-rays;
- panoramic x-rays;
- cephalometric x-rays ;
- intraoral video camera images;
- other types of images.

It is also possible to choose to save the treated or original image data, with or without burnt-in annotations.

If treated image data is selected then the image is stored as displayed in VixWin.

If annotations is selected, then any annotations applied to the image, will be burnt into the image data before storing the image. Applying annotations to images is described in the VixWin Platinum User Manual.

Saving the Configuration changes

Click **SAVE** to save the configuration.

Click **EXIT** to exit the configuration window.

NOTE: If the VixWin program is running, then it must be restarted before configuration changes take effect.

CHAPTER 3 Operation with Remote Servers


Modality Worklist Server

A Modality Worklist Server for a RIS (Radiology Information System) can be interfaced with VixWin to generate a list of patients/studies.

NOTE: VxDICOM must be configured to use the Worklist server by setting the Patient Database to DICOM Worklist.

The **VixWin DICOM Worklist SCU** window can be displayed at any time by clicking on the patient icon button, or by choosing **File > Patient Worklist**.

When the VixWin DICOM Worklist SCU window is open, the user can then click one of the buttons below to query for the patient list.

Worklist: Query for the complete Worklist.

Find by patient: Query the Worklist by Patient, Patient ID or Accession Number.

Find by date: Query the Worklist by Date.

NOTE: The query will also include the Scheduled Station Name DICOM tag (0040,0010). This tag can be configured in VixWin by choosing **Options > Preferences > Customize** and entering the Station Name.

Double click a patient from the list and the selected patient name appears in the VixWin program title bar.

When passing from one patient to another VixWin closes all displayed images. This ensures that all displayed images belong to the currently selected patient.

We recommend starting a work session by selecting a patient from the Worklist, before acquiring images.

Image Storage

A DICOM Storage Server of a PACS (Picture Archiving and Communication System) can be interfaced with VxDICOM for storing images acquired in VixWin.

Images can be saved to PACS by clicking the **Send to PACS** button or by choosing **File > Save to DICOM Server**.

NOTE: If you have set the VixWin Image Database to DICOM then database images are saved to PACS as default. Images are saved automatically in Full Mouth view, while in Windows Layout view you may also use the **File > Save Image** command or button.

Before saving an image to PACS, you can specify DICOM image data to be saved with the image by choosing **File > DICOM Dataset**.

See the VixWin Platinum User Manual for assigning tooth numbers, notes etc to images. Select the Image tab to set the Positioner Type, Patient Orientation and Image Laterality. Remember to assign data before sending the image to the DICOM Server. Once the image is stored it cannot be changed.

Query/Retrieve

A DICOM Storage Server (PACS) can be interfaced to VxDICOM for retrieving previously stored images.


Images can be retrieved from PAC by clicking the PACS Query/Retrieve button or by choosing **File > DICOM Query/Retrieve**.

NOTE: If you have set the VixWin Image Database to DICOM then database images are retrieved from PACS as default and so you can also use the **File > Load Image** command or button.

If a patient is selected in VixWin the Q/R window will query for this patient.

To query the PACS click the Find button in the Q/R module window to display the search dialog. Type in the Patient Name, Patient ID and/or other search parameters and click the Search button to start the query.

The list of patients, studies and series/images that match the search is displayed.

You can search through the patient folders to view the studies with corresponding series/images. Patients are displayed together with name, ID, date of birth and sex. Studies are displayed with date and ID. Images/ Series are shown with corresponding image type.

In order to retrieve images, right-click the patient, study or series/image and choose Retrieve.

The different image types are depicted by icons, as follows:

IO image: Intraoral x-ray images.

PX image: Panoramic x-ray images.

DX image: Cephalometric x-ray images

GM/ES/XC images: Dental Color camera images.

Generic image

Patient Information

Patient Information

Patient ID: 1234ABCD

Last name: Doe

First name: John

Date of birth: 14 MAY 1961

Sex: ☒ Male ☐ Female

OK Cancel

Patient Information can be manually entered by choosing **File > Patient Information**.

This is often useful if a Worklist is not available for selecting the patient.

DICOM Removable Media

The VixWin Platinum software also includes DICOM Removable Media functionality through the DICOM Import and DICOM Export menus. This functionality is included in the standard VixWin Platinum software and is therefore independent of the VxDICOM module. This manual describes the DICOM networking functionality only. The DICOM Import/Export functions are described in the VixWin Platinum User Manual


DICOM Print

Images can be printed from VixWin DICOM to DICOM printers.

NOTE: VixWin DICOM only prints grayscale images (x-rays). Color images are not supported.

- Choose **File > Print Images DICOM** to print the displayed images in VixWin Platinum Windows Layout viewing mode.
- Choose **File > Print Image DICOM** to print the selected displayed image in VixWin Platinum Full Mouth viewing mode.
- Choose **File > Print Mount DICOM** to print the displayed mount in VixWin Platinum Full Mouth viewing mode.

NOTE: The individual images in Print Mount DICOM will not be scaled to real dimensions.

Select an image and then choose **File > Print Images DICOM** to open the **Print DICOM** dialog box above.

The various settings are described below (“*” indicates the default values):

1. Priority
 - HIGH *
 - MED
 - LOW
2. Medium Type
 - PAPER
 - CLEAR FILM
 - BLUE FILM *
3. Film Destination
 - PROCESSOR *
 - MAGAZINE
4. Annotations (DICOM Information)

VixWin DICOM can print DICOM information with the image in the form of annotation text:

 - NONE * (Do not print annotations)
 - DICOM (Print annotations)
 - VIXWIN (Obsolete. Backward compatibility only)

Consult the DICOM printer Conformance Statement for the configuration parameters of your printer model. Typical parameters are:

DICOM Printer	Display Format ID	Position
AGFA	ANNOTATION	1 or 6
FUJI	FORMAT1	1 or 6
KODAK	LABEL	1
CODONICS	1T	1

5. Film Format

Specifies how images are positioned on the film. Use 1 X 1 for printing single images on a page. Default: 1 X 1
6. Film Size

The size of the film to be used. Default: 8IN X 10IN
7. Orientation

LANDSCAPE or PORTRAIT *

8. Trim Around Image
Should a trim be printed around each image on the film:
YES
NO *
9. Border Density
Border color of the areas around the image boxes on the film.
BLACK *
WHITE
10. Empty Image Density
Density of the image box area on the film, not containing the image.
BLACK
WHITE *
11. Min Density
Minimum density of the images on the film, expressed in hundredths of OD. Default value: 0.1
12. Max Density
Maximum density of the images on the film, expressed in hundredths of OD. Default value: 3.5

CHAPTER 4

Glossary

Below is listed a glossary of terms used in this manual.

AE (Application Entity) - a DICOM conformant Application on the network.

Calling AE Title - the name (title) of the SCU that is attempting to start an Association with a SCP. This is the name of the DICOM Application that does the calling.

Called AE Title - the name (title) of the SCP with which a SCU attempts to start an association. This is the name of the DICOM Application that is being called.

Conformance Statement - A document whose organization and content are mandated by the DICOM Standard, which allows users to communicate how they have chosen to comply with the Standard in their implementations.

DICOM - Digital Imaging and Communications in Medicine.

DX - Digital Radiology

ES - Camera Endoscopic

GM - General Microscopy

HIS - Hospital Information System

Host Name - This is a friendly name for the computer. It can be used in place of the IP address if the computer is local. If the computer is accessed over the internet, its IP address needs to be used

IO - Intraoral Radiography

IP Address - A unique identifier for the network interface of a computer on a TCP/IP network. An IP address is typically comprised of four octets, separated by dots (.),

with each octet capable of representing a number from 0 to 255. For example:
192.168.10.1

PACS - Picture Archiving and Communication System

Port - Port address used for communication between the SCP and SCU computers.

PX - Panoramic X-Ray

RIS - Radiology Information System

SCP (Service Class Provider) - A device which provides DICOM services which are utilized by another device (SCU).

SCU (Service Class User) - A device which utilizes the DICOM Services which are provided by another device (SCP).

Modality Worklist Server - A device used to download worklists, containing patient demographic data, from hospital information systems (HIS/RIS) to the acquisition modality (VixWin DICOM)

XC - External Camera Photography

Applicable Standards

IEC 601.1.4 (Medical Electrical Equipment. Collateral standard: Programmable electrical medical systems).

93/42 EEC, European Directive concerning medical devices

Digital Imaging and Communication in Medicine (DICOM) standard, V3.0

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Kaltenbach & Voigt GmbH