VISIODENT

V I S I O C A M 3

USER MANUAL
MANUEL D'UTILISATION

Contents

| | | |
|-----------|---|----|
| 1 | FOREWORD | 2 |
| 2 | DENTAL CAMERA INTRODUCTION | 3 |
| 3 | SAFETY INSTRUCTIONS | 4 |
| 4 | REGULATORY REQUIREMENTS | 5 |
| 5 | INSTALLATION | 10 |
| 6 | CONNECTING TO A COMPUTER | 11 |
| 7 | DESCRIPTION OF THE CONNECTION BOXES | 12 |
| 8 | MAINTENANCE | 14 |
| 9 | AFTER-SALES SERVICE | 16 |
| 10 | TECHNICAL FEATURES | 19 |

FOREWORD

Thanks for trusting us by purchasing this device.

In order to optimize the use of the device, while taking all the necessary precautions, it is recommended to carefully read and follow the user manual.

Always carefully consider the messages CAUTION, WARNING, and NOTE when using the system.

CAUTION: *the term CAUTION describes potential incidents likely to jeopardize safety.*

WARNING: *the term WARNING refers to the incidents likely to disturb the smooth running of the imaging system.*

NOTE: *the term NOTE highlights particular points in order to facilitate the system maintenance or to clarify important information.*

DENTAL CAMERA INTRODUCTION

This is an intra-oral camera designed for dental applications. It allows visualization of anatomical and pathological details that cannot be seen with the naked eye, as well as the pathology control before and after treatment. Its extra-light handpiece, sensitivity and resolution performance, with respect to shades, make it an ideal medical tool.

This dental camera is composed of a camera handpiece, a connection box as well as various accessories necessary for it to work.

VISIOCAM 3

- 10 hygienic protection covers.
- 1 handpiece integrating the camera electronics and lighting.
- 1 handpiece holder.

DOCK USB2

- A connection box with a 3.5 metre cable to link the handpiece to the computer USB port.
- Imaging software: VISIODENT Imaging.
- An user's manual included in the Imaging software.
- A quick start guide.

This equipment was packaged in a custom carton. This carton should be kept for possible transportation. As a complement of the dental camera, some dental barriers necessary to the intra-oral use of the camera are included.

For more details about these products, please refer to our catalog or contact our commercial service.

NOTE:

The device was designed and developed with its accessories in order to guarantee to you safety and maximum performance. The use of accessories originating from any other company than Visiodent can represent a risk for you, your patients and your device.

SAFETY INSTRUCTIONS

- Do not expose the dental camera to water spray and do not store it in a humid environment (electrocution risk).
- Install the dental camera in a clean, dry, and well-ventilated place.
- Disconnect the camera from the power supply if you are not going to use it for several days. Do not pull on the cable.
- Never compress or nip the handpiece cable.
- Never expose the product to high vibrations.
- Do not drop the camera.
- This camera should in no way be immersed in any liquid, nor should it be autoclaved.
- The surface temperature in the light emission area can be just above 41°C (after several minutes using it). So, avoid maintaining this emission area in contact with the patient's mouth.
- The VISIOCAM 3 camera is a product using class II LEDS according to IEC 60825. Do not stare at them in order to avoid any ocular risk.
- For each new patient, it is essential to use the hygienic protection covers provided with the camera or provided as a complement.

**NOTE:**

If the hygienic protection is torn while examining a patient or if the handpiece was "infected" while withdrawing the hygienic protector, it is essential to totally disinfect the handpiece. In order to do this: please refer to the maintenance chapter.

REGULATORY REQUIREMENTS

4.1. COMPLIANCE WITH STANDARDS AND REGULATIONS

This product was designed and manufactured by a company having an authorized quality foundation. It meets the requirements of the European directive 93/42/EEC, relative to medical devices. Therefore, it particularly meets electrical safety and electromagnetic compatibility standards (IEC) (CEM).

4.2. ELECTROMAGNETIC INTERFERENCE AND ELECTROSTATIC DISCHARGES

Electromagnetic compatibility (CEM) is the ability of electronic device elements to correctly interact in an electronic environment. Although this dental camera system was designed according to this compatibility and complies with the electromagnetic interference thresholds established by the regulatory agency, there is no guarantee about the interference likely to occur on a particular installation.

If the device generates interference with radio communication services (which can be determined by switching it off and on), the user is recommended to try to correct this phenomenon by taking in whole or part in of the following measures:

- Change the receiving antenna orientation.
- Reposition the product according to the receiver.
- Move the receiver further away from the computer.

VISIOCAM 3 dental camera system is designed and tested to be used in a home environment, class B Group 1, according to CISPR11 standard.

4.3. MEDICAL DEVICE VIGILANCE

As any medical device, this device is subjected to the medical device vigilance dispositions; any serious dysfunction should then be the subject of a description to the competent authorities and to the manufacturer as soon as possible and as precisely as possible.

4.4. END OF LIFE

This device bears the recycling symbol according to the European directive 2002/96/EC for electric and electronic equipment waste (DEEE or WEEE). Properly disposing of this device will contribute to avoiding any further damage to the environment or human health.

The symbol  present on the device or on the accompanying documentation indicates this product cannot be in any case treated as household waste. Therefore, it should be given to a waste collection centre in charge of electric and electronic equipment recycling.

For scrapping, respect the standards relative to waste disposal in force in the installation country. For more details about the device treatment, recuperation and recycling, please contact your dental device distributor (or failing that, the VISIODENT website), so that you can be provided the proper procedure.

4.5. ELECTROMAGNETIC COMPATIBILITY

| Guide and declaration of the manufacturer - electromagnetic emissions | | |
|---|---------------|--|
| VISIOCAM 3 device is intended to be used in the electromagnetic environment specified below. The user should make sure it is used in this environment. | | |
| Emission trial | Compliance | Electromagnetic environment - Guide |
| RF emissions CISPR 11 | Group 1 | VISIOCAM 3 device only uses radio energy for its internal functions. Therefore, its RF emissions are very low and are not likely to cause interference with electronic devices nearby. |
| RF emissions CISPR 11 | Class B | |
| Harmonic emissions IEC 61000-3-2 | No applicable | |
| Voltage fluctuations / Flicker IEC 61000-3-3 | No applicable | |

| Guide and declaration of the manufacturer - electromagnetic immunity | | | |
|---|--|--|---|
| VISIOCAM 3 device is intended to be used in the electromagnetic environment specified below. The user should make sure it is used in this environment. | | | |
| Immunity trial | CEI 60601 Severity level | Compliance level | Electromagnetic environment Guide |
| Electrostatic discharges IEC 61000-4-2 | ± 6 kV when in contact ± 8 kV in the air | ± 6 kV ± 8 kV | The floor should be wooden, concrete or tile. If the floor is covered with a synthetic material, the relative humidity should be at least 30%. |
| Far transient bursts IEC 61000-4-4 | ± 2 kV for the feed cables ± 1 kV for the input/output cables | ± 2 kV ± 1 kV | The main power supply quality should be one of a traditional commercial or hospital environment. |
| Voltage shocks IEC 61000-4-5 | Differential mode ± 1 kV Common mode ± 2 kV | ± 1 kV N.A. | The main power supply quality should be one of a traditional commercial or hospital environment. |
| Dips, brief outages and power voltage variation IEC 61000-4-11 | <ul style="list-style-type: none"> $<5\% U_T$ - for 10 ms $40\% U_T$ - for 100 ms $70\% U_T$ - for 500 ms $<5\% U_T$ - for 5 s | $<5\% U_T$ 10 ms $<40\% U_T$ 100 ms $<70\% U_T$ 500 ms $<5\% U_T$ 5 s | The main power supply quality should be one of a traditional commercial or hospital environment. If the user of VISIOCAM 3 device requires it to continue to operate during main power supply outages, it is recommended VISIOCAM 3 device is fed by an inverter or a battery. |
| Magnetic field with the network frequency (50/60 Hz) | 3 A/m | 3 A/m | The magnetic field with the network frequency should be at a characteristic level of a location in a traditional commercial or hospital environment. |
| Note: U_T is the power voltage nominal value applied during the trial. | | | |

| Guide and declaration of the manufacturer - electromagnetic immunity | | | |
|--|--------------------------------|---------------------|---|
| VISIOCAM 3 device is intended to be used in the electromagnetic environment specified below. The user should make sure it is used in this environment. | | | |
| Immunity trial | CEI 60601 Severity level | Compliance level | Electromagnetic environment Guide |
| Conducted RF IEC 61000-4-6 | 3 Vrms 150 kHz to 80 MHz | 3V | Portable and mobile RF communication devices should not be used at a distance from VISIOCAM 3 device, including the cables, lower than the recommended separation distance, calculated with the applicable formulas depending on the emitter frequency. Recommended separation distance $d = 1,16\sqrt{P}$ |
| Radiated RF IEC 61000-4-3 | 3 V/m 80 MHz to 2.5 GHz | 3V/m | $d = 1,16\sqrt{P}$ 80 MHz to 800 MHz $d = 2,33\sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum rated output of the transmitter in watts (W) by the transmitter manufacturer and d the recommended separation distance in metres (m). The field levels emitted by the fixed RF transmitters, determined by an electromagnetic measurement of the site ^a , should be lower than the compliance level in each frequency band. Interference may occur in the vicinity of the devices bearing the following symbol:  |
| Note 1: At 80 MHz and 800 MHz, the higher frequency band applies. Note 2: These recommendations may not apply in every situation. Electromagnetic wave propagation is modified by the absorption and reflection due to the structures, objects and persons. | | | |

- a The fixed transmitter field levels, such as the base stations of the radio telephones (cellular/wireless) and the terrestrial mobile radios, amateur radio, AM, FM, and TV radio communication cannot be theoretically assessed precisely. To obtain the electromagnetic environment due to the fixed RF transmitters, a site measurement should be performed. If a field level measured in the use environment of VISIOCAM 3 device exceeds the compliance levels above applicable, the good operation of VISIOCAM 3 device should be checked. If abnormal operations are proved, some further measures should be taken, such as reorientation or relocation of the standard device.
- b Above the 150 kHz to 80 MHz frequency band, the field level should be lower than 3 V/m.

Recommended separation distances between the portable and mobile RF communication devices and VISIOCAM 3 device

VISIOCAM 3 device is intended to be used in an electromagnetic environment in which the irradiated RF disturbances are checked. The user of VISIOCAM 3 device can help to avoid electromagnetic interference by maintaining a minimal distance between the portable and mobile RF communication devices (transmitters) and the recommended VISIOCAM 3 device such as recommended below, depending on the maximal output power of the communication device.

| Rated maximal output power of the transmitter W | Separation distance depending on the transmitter frequency m | | |
|--|---|--------------------|--------------------|
| | 150 kHz to 80 MHz | 80 MHz to 800 MHz | 800 MHz to 2,5 GHz |
| | $d = 1,16\sqrt{P}$ | $d = 1,16\sqrt{P}$ | $d = 2,33\sqrt{P}$ |
| 0.01 | 0.116 | 0.116 | 0.233 |
| 0.1 | 0.366 | 0.366 | 0.736 |
| 1 | 1.16 | 1.16 | 2.33 |
| 10 | 3.66 | 3.66 | 7.36 |
| 100 | 11.6 | 11.6 | 23.3 |

For the transmitters whose maximal output is not listed above, the recommended separation distance d in metres (m) can be determined by using the equation applicable to the transmitter frequency, where P is the maximal output of the transmitter in watts (W) rated by the transmitter manufacturer.

Note 1: At 80 MHz and at 800 MHz, the separation distance given in the higher frequency band applies.

Note 2: These recommendations may not apply in every situation. The electromagnetic wave propagation is modified by absorption and reflection due to the structures, objects and persons.

INSTALLATION

5.1. COMMISSIONING

Fastening the handpiece holder:

1. Choose an area without obstructions that can be easily accessible for use;
2. Use the wipe provided to clean the surface on which you are going to fasten the holder;
3. Remove the double-sided tape adhesive protection that is on the holder, place it and press it several times to correctly fix. The maximum sticking performance is obtained after two hours so avoid any stress on the holder during this period.

CAUTION:

This holder is equipped with magnets that can damage devices sensitive to magnetic fields. Make sure you do not install this holder near these devices (cathode ray tube video screen, magnetic videotapes, etc.).

5.2. FURTHER CONNECTION BOXES (optional)

You can install a connection box near each dental chair (no limitation). You will just have to transport the handpiece from one chair to the other.

The handpiece holder is intended to maintain the connecting cable connector when the cable is not linked to the handpiece.

When you disconnect the handpiece connecting cable to take it to another chair, or when you put the handpiece on its holder, the last stored image displays on the screen* (or the last four if you were in four-image mode, or the color-bar pattern if no image was stored).

*Except on Dock USB2.

CONNECTING TO A COMPUTER

6.1. REQUIRED CONFIGURATION FOR THE COMPUTER

To use the VISIOCAM 3 system, you must make sure the computer and its peripherals do not have any use limitation that could concern life safety. It should also meet the following requirements:

| | Minimal configuration | Recommended configuration |
|-------------------|---|--|
| Operating system | Windows 2000 Pro SP4 - XP Pro SP1 | Windows XP Pro SP3 - VISTA Business SP1 |
| Processor | Intel® Pentium III - 500 MHz | Intel® Core 2 |
| Memory | 512 MB | 2 GB |
| Hard disk | 20 GB | 160 GB |
| USB ports | 2 x USB2.0 Hi-Speed ports | 4 x UB 2.0 Hi-Speed ports |
| Video board | Graphic board 32 MB of unshared video RAM compatible with DirectX 9 | Graphic board with Nvidia chipset or exclusive video ATI / 256 MB of unshared video RAM compatible with DirectX 9 or higher. |
| USB Chipset | Intel or NEC | Intel or NEC |
| Screen resolution | 1024 x 768 | 1024 x 768 or higher |
| Standard | IEC 60950 | IEC 60601-1 |

6.2. DOCK USB2 CONNECTION

- Connect the USB cable between the connection box and one of the computer USB ports.
- Connect the connecting cable to the handpiece.

6.3. VISIODENT SOFTWARE CONFIGURATION

- Insert the VISIODENT CD-ROM into the driver. It automatically executes.

NOTE :

Refer to the installation manual of VISIODENT software that is on the CD-ROM VISIODENT Imaging.

6.4. VISIODENT SOFTWARE CONFIGURATION

Refer to the installation and configuration manual of the VISIOCAM 3 camera that is on the VISIODENT CD-ROM in the document directory.

DESCRIPTION OF THE CONNECTION BOXES

7.1. DOCK USB2 POWER SUPPLY

The dental camera electrical supply is directly performed through the computer USB port. The voltage powering the camera is of continuous 5 V low voltage type  (0.5 A).

7.2. USB OUTPUTS

This camera has a digital USB 2.0 output that can be connected to the USB port2 of a computer.

7.3. IDENTIFICATION

The indications born on the boxes identify the dental camera according to the international standards CEI 601-1 and CEI 417.

 Continuous voltage.

USB  USB2 output.

“BF type camera”.

Caution, please consult the ACCOMPANYING DOCUMENTS.

Disposal of electric and electronic equipment marketed after 13/August/2005. This symbol indicates that the product cannot be treated with domestic waste

For active medical devices, this symbol is associated to the manufacturing year (expressed with four digits).

For active medical devices, this symbol is associated to the manufacturer name and address.

Product compliance according to the European directive 93/42/EEC relative to medical devices.

Grounding.

The devices that connect to video or USB outputs should comply with IEC 60950 standard.

MAINTENANCE

For each new patient, it is essential to use the hygienic protection covers provided with the camera or provided as a complement.

Before using the system for the first time, it is imperative to follow the complete disinfecting procedure.

Any intra-oral camera returning from service or maintenance should also receive the complete disinfecting procedure before being used.

! WARNING :

Do not use products containing:

- *Ammoniac, trichloroethylene*
- *Dichloroethylene*
- *Ammonium hydrochlorid*
- *Chlorinated and aromatic hydrocarbon*
- *Ethylene dichloride*
- *Methylene chloride*
- *Ketones*

Deterioration risks of the plastic parts.

! WARNING :

Do not directly spray disinfecting products on VISIOCAM 3 intra oral camera.

8.1. HANDPIECE AND CONNECTION BOX MAINTENANCE:

| DESCRIPTION | RECOMMENDATIONS | USE INSTRUCTIONS AND PRECAUTIONS | |
|---------------|---|--|---|
| | | ✓ | ✗ |
| Disinfecting. | <ul style="list-style-type: none"> • Surface cleaning and disinfecting wipes of Wipes Septol™ type - Pierre Rolland. | <ul style="list-style-type: none"> ✓Take the cleaning cloth, wring it, and then swipe the equipment until obtaining visible cleanliness. ✓Allow to dry in the open air. ✓Carefully close the packaging box. | <ul style="list-style-type: none"> ✗ Do not scrub ✗ Do not rinse. ✗ Do not immerse in a disinfecting liquid. |

NOTE :

In case of contact with blood or important soiling, it is strongly recommended to follow a disinfecting process. First of all, clean the handpiece with disinfecting wipes, then roll up several disinfecting wipes around the handpiece and leave on for 15 minutes.

AFTER-SALES SERVICE

9.1. WARRANTIES

VISIODENT ensures its products to be free from material and manufacturing defects for a period of one (1) year beginning on the purchase date. This warranty does not apply to misused, modified, untended, or accidentally damaged products, or products subject to abnormal use and handling conditions. The distributors are not authorized to extend any warranty on behalf of VISIODENT.

The entire liability of VISIODENT is limited to its convenience when replacing or repairing, free of charge, the defective product if it has been sent to VISIODENT After-Sales Service during the warranty period only.

Outside of France, the access to the warranty is only possible if the product was bought from one of the points of sale authorized by VISIODENT in the country where it will be used.

THIS WARRANTY CONSTITUTES THE ONLY AND UNIQUE REMEDY. IT REPLACES ANY OTHER WARRANTY, FOR EXAMPLE, OF ADEQUACY TO A PARTICULAR AIM, SHOULD IT BE EXPLICIT OR IMPLICIT. VISIODENT SHALL NOT BE LIABLE FOR ANY PARTICULAR DAMAGE, INDIRECT, ACCIDENTAL OR CONSEQUENTIAL, OR FOR ANY DETERIORATION OR DATA LOSS, ON A CONTRACTUAL, NON-CONTRACTUAL OR OTHER BASIS.

The liability exclusion or limitation for direct or indirect damage does not apply under the regulatory or legal rules in force in some countries and the present exclusion may not apply to a purchaser in those countries.

9.2. IN CASE OF FAILURE

| PROBLEMS | CAUSES | SOLUTIONS |
|---|--|---|
| No image displays on the screen and the camera LEDs are not on. | <ul style="list-style-type: none"> • Defective power supply. • Connection problem. | <ol style="list-style-type: none"> 1. Check the power supply is correctly connected to the network and to the connection box. 2. Check the connecting cable is correctly connected to the handpiece and to the connection box. |
| The camera switches on, but no image displays on the screen. | <ul style="list-style-type: none"> • Configuration. • Driver. • Connection problem. | <ol style="list-style-type: none"> 1. Check the camera is correctly set up in VISIODENT (please, refer to VISIODENT owner's manual). 2. Check the camera is correctly detected in the device driver (correct installation of its driver). 3. Check the USB cable coming from the DOCK is correctly connected to the HUB. |
| An image displays on the screen, but the quality is not satisfactory. | <ul style="list-style-type: none"> • Camera driver configuration. | Check the camera configuration in the VISIODENT software (brightness, contrast, saturation, etc.). Please refer to VISIODENT owner's manual. |
| An image displays, but it is not very clear (blurry) | <ul style="list-style-type: none"> • Hygienic protector. | Check the hygienic protector is correctly positioned on the camera head. |

If the defect persists and you will return the dental camera to your distributor, be sure to send it in its original packaging. In the same way, the dental camera should be sent to your distributor in its totality (Connection box, handpiece, cables). Then, please, enclose with your packing list, a brief note relative to the noticed defect.

If some parts constituting the camera happened to break, it is imperative to include everything so that the defective parts can be replaced.

When your material is returned to you, you should check its condition and note any reservation on the delivery slip, if necessary. You then have 48 hours to confirm by registered letter sent to the carrier. After this period, the carrier will be able to deny these reservations.

If any material returned to you is damaged during transportation, the repair charges will be billed to either the carrier (if the reservations were made within the period) or the recipient. So confirm as soon as possible that all materials are working properly.

TECHNICAL FEATURES

VISIOCAM 3

- High sensitivity CCD 1/3".
- Resolution: (752 x 582) PAL; (768 x 494) NTSC.
- Definition: 470 lines.
- Sensitivity: 2 lux.
- Ratio Signal/Noise: 52 dB.
- Electronic diaphragm: automatic (1/50 to 1/100 000).
- Lighting: 8 LED.
- Freeze image  or foot control (optional)
- Adjustments: automatic.
- Non-inverted image.
- Visual field: 108°.
- Handpiece dimensions: L: 198; W: 27; H: 18 mm.
- Distal part dimensions: W: 16.5 x H: 12 mm.
- Handpiece weight: 50 g.
- applied part of BF type.

DOCK USB2

- Cable length: 3.5 m.
- 1 USB output 2.0.
- Dock dimension: L: 100; W: 46; H: 20 mm.
- Dock weight: 165 g.
- Operating temperature: +10 °C to +40°C.
- Storage temperature: -20°C to +45°C.
- Relative humidity: 10 % to 90 %.
- Atmospheric pressure: 900 hPa to 1,060 hPa.
- Continuous service.
- Not protected against water chutes (IPX0).
- Not adapted to the use in presence of an anaesthetic mixture flammable with air, oxygen or dinitrogen monoxide.
- Complies with the European directive 93/42/EEC.
- Complies with EC60601-1 standard.
- Complies with IEC60601-2-18 standard.
- Complies with UL 60601-1 et CSA 60601-1 standard.

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