MEMORY *m*

VIEWWIRELESS

WIRELESS

USB

RSV-HD

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
CONTENTS

1 SAFETY INSTRUCTIONS AND REGULATIONS

1 1 System maintenance	49
1 2 Compliance with standards	49
1 3 Workstation compliance	49
1 4 Marking on the controller box	50
1 5 Operational limitations	50
1 6 List of the elements contained in your system	50
1 7 Minimal configuration	51

2 RSV-HD VIEWWIRELESS

2 1 Installation and configuration	52
2 2 Use	59
2 2 1 Wireless connection with the computer	59
2 2 2 Use in autonomous mode	61
2 3 Batteries	63
2 4 Caution in the use of your sensor	64
2 5 Caution in the use of your controller box	65
2 6 Technical characteristics	66
2 7 Fixation	66
2 8 Information concerning the RSV-HD Viewwireless	67

3 RSV-HD MEMORY

3 1 Installation and configuration	68
3 2 Use	70
3 2 1 Wireless connection mode with the computer	70
3 2 2 Use in autonomous mode	72
3 2 2 1 Classic stand alone mode	72
3 2 2 2 Advanced stand alone mode	74
3 2 2 3 Transfer mode	75
3 3 Batteries	77
3 4 Caution in the use of your sensor	77
3 5 Caution in the use of your controller box	77
3 6 Technical characteristics	77
3 7 Fixation	77
3 8 Information concerning the RSV-HD Memory	77

4 RSV-HD WIRELESS

4 1 Installation and configuration	78
4 2 Use	78
4 3 Batteries	79
4 4 Caution in the use of your sensor	79
4 5 Caution in the use of your controller box	79
4 6 Technical characteristics	79
4 7 Fixation	80
4 8 Information concerning the RSV-HD Wireless	80

5 RSV-HD USB

5 1 Installation	80
5 2 Configuration	81
5 3 Use	82
5 4 Caution in the use of your sensor	82
5 5 Caution in the use of your controller box	82
5 6 Technical characteristics	82

6 TECHNICAL CHARACTERISTICS OF THE CCD SENSOR

7 KIT OF POSITIONNERS

8 WARRANTY

WELCOME

You have recently purchased a digital radiology system from Visident group. We thank you for your trust and we will do our best to make you satisfy.

This user manual describes the way to use the new RSV-HD line.

Visident company shall not be held responsible for the consequences resulting from non respect of all the precautions described in the present manual.

Visident, Visident XP, RSV - Radiology System Visident -, RSV-HD, Visident Imaging, as well as logos are Trademarks of the Visident Group. Microsoft, Windows are Trademarks of the Microsoft Corporation.

All the other product names mentioned in the present manual are used for identification purposes and may be brands or trademarks registered by their owners.

This manual is the intellectual property of the Visident Group with all rights reserved. Any partial or complete reproduction of this manual in any form is strictly prohibited without the prior written consent of the Visident Group. Because of the constant evolution of its products, Visident reserves the right to modify the contents of this manual even without prior notice.

The RSV-HD technologies are covered by international patents by Visident Group.

1 SAFETY INSTRUCTIONS AND REGULATIONS

1 | 1 System maintenance

Your digital radiology system Visident must be installed and used respecting all the necessary hygienic and safety instructions (for both practitioner and patient) contained in this manual.

Your digital radiology system Visident is only to be used only in areas or rooms where regulations regarding electrical safety in medical premises are strictly enforced.

This device has to be installed in a dry place free of humidity. Any liquid penetration in the controller box or in the sensor could cause damage, provoke a short circuit or corrosion and prevent the system from working properly.

Your digital radiology system Visident should always be disconnected before maintenance work or disinfection is carried out.

As they evaporate, some disinfection products form explosive or inflammable mixtures. If this type of disinfectant is used, wait for its complete dispersion before using the system again.

Finally, it is advised to transport the system in its original package.

1 | 2 Compliance with standards

The RSV-HD line complies with EC Medical Devices Directive 93/42/EEC on medical instruments. According to the degree of harmful liquid penetration, it is considered to be a standard instrument.


The line complies with European standard applying to this case: general Rules of Safety and Electromagnetic Compatibility (NF EN IEC 60601-1 and NF EN IEC 60601-1-2) as well as with the medical safety and electromagnetic standard IEC 601-1.

In accordance with the obligations provided for the EC label, the user is bound to declare to the Ministry of Health or the proper Health Authority any relevant information on incidents related to the device as well as any alteration of the specifications of the device that could cause the death or the worsening of the state of health of the patient or of the person who manipulates the device.

At the same time, Visident shall be kept informed of this declaration so that it can perform the obligation provided for the Directive.

1 | 3 Workstation compliance

Before installing your digital radiology system Visident, please read and apply the instructions shown on the diagram carefully.

The computer and display monitor (non-medical electrical equipment) **must be situated outside the area shown here** (1,5 meters between the patient and parts of the system which can unintentionally be touched or between patient and other persons touching parts of that system).

The computer with which the digital radiology system Visident is installed, and all additional devices connected; must have the CE marking (IEC/EN 60950). If not, the equipment must be connected to an isolating transformer (IEC 60989 - withstand 1500 V).

In order to have the best image quality, it is advised to position the screen in order to avoid direct light reflections (interior or exterior light).

1 SAFETY INSTRUCTIONS AND REGULATIONS

1 | 4 Marking on the controller box

This label guarantees the compliance of the system with EC Medical Directive 93/42/EEC on medical instruments.

Read the enclosed instructions (in this Manual).

This device bears the recycling symbol according to the European directive 2002/96/EC for electric and electronic equipment waste (DEEE or WEEE). Properly disposing of this device will contribute to avoiding any further damage to the environment or human health.

Non-ionizing radiation

SN Serial number

Manufacturer.

1 | 5 Operational limitations

The device is designed for operation within the following limits :

Non-operating environment :

- Temperature range: -40 to + 70°C
- Humidity range : 10 to 95% RH

Operating environment :

- Temperature range : 10 to 40°C
- Humidity range : 20 to 60% RH
- Liquid penetration : index of penetration IPX0
- Atmospheric pressure: no notable influence in normal operating conditions. 500 hPa to 1060 hPa during transportation and storage.

Do not touch the computer screen and the sensor simultaneously ; this could result in severe damage to the sensor.

1 | 6 List of the elements contained in your system

For RSV-HD Memory, Viewwireless and Wireless :

- 1 Controller box
- 2 kinds of fixations : 1 scratch strap fixation and 2 plastic straps fixation
- 1 fixation support
- 8 rechargeable batteries AAA (4x2)
- 1 batteries charger to connect directly to the system
- 1 external batteries charger
- 4 adapters for US, europeans (UK - Europe) and Australian plugs
- 1 Sensor Size 1 or 2
- 1 Positionners kit Size 1 or 2.
- 1 USB Key Bluetooth
- 1 CD-ROM containing :
 - VisiDent Imaging software
 - VisiDent Imaging user manual
 - 1 VisiDent Imaging installation manual
 - This RSV-HD user and installation manual

For RSV-HD USB :

- 1 Controller box
- 2 kinds of fixations : 1 scratch strap fixation and 2 plastic straps fixations
- 1 fixation support
- 4.8 meters USB Cable
- 4 adaptater for US, europeans (UK - Europe) and Australian plugs
- 1 sensor Size 1 or Size 2
- 1 Positionners kit Size 1 or 2
- 1 CD-ROM containing :
 - VisiDent Imaging software
 - VisiDent Imaging user manual
 - 1 VisiDent Imaging installation manual

1 | 7 Minimal configuration

The computer

The computer plays a vital role since it is responsible for storing the digital X-rays (which are heavy files requiring a considerable amount of disk space) and running the Visident Imaging software.

Its characteristics :

Intel Processor Pentium 3.1 GHz Minimum
256 Mb Ram
Hard disk 40 Gb Minimum
CD-ROM player
USB 2 port compatible USB1 (for RSV-HD)
Video card AGP type 32Mb RAM
OS Windows XP Professional, Mac OSX

The computer must comply with standard IEC/EN 60950.

The Display Monitor

Its role is also very important because it allows you to visualise the digital X-rays in an optimal way.

It is recommended to use multisync colour monitor with a minimum display resolution of 1024 x 768 (this resolution can go up to 1280 x 1024).

It is also recommended to use a 17inch monitor with a dot-pitch inferior or equal to 0.26 microns.

The X-ray tube

The RSV-HD System is compatible with any kind of x-ray tube if the x-ray tube is able to provide an necessary x-ray dose and time exposure.

However, it is recommended to use x-ray tube with a voltage between 50 and 70 kV and an intensity from 6 to 8mA. The timer should be able to offer an exposure time going from 60 to 350ms. For any other information please refer to the Visident Imaging user manual.

② RSV-HD VIEWWIRELESS

2 | 1 Installation and configuration

Visident Imaging installation


To use the RSV-HD system you must install the latest version of Visident Imaging. For this, use the installation CD provided with the system and follow the instructions. your RSV-HD System is provided with its own « Visident Imaging software ». However, it is also possible to use this system with other imaging softwares thanks to some bridging tools that can be provided by Visident upon demand. Visident takes no responsibility for the operation of software other than Visident Imaging.

Installation of the Bluetooth key on the computer

This key enables your RSV-HD system to communicate with the computer. You will find the way to install the key :


Insert the Bluetooth key provided in the USB port of your computer.

A window appears to install the new peripheral detected, but we will not use this assistant.


To install this peripheral, you will use the driver provided in the Visident Imaging CD. Browse this CD and look for the file « Driver Bluetooth USB Adapter ». Open this file and double click on « Setup.exe ».


Then, the installation program of the Bluetooth key starts. The following window appears :


Click on « Next> ». A new window appears :


Please, read the end user license agreement and select « I accept the terms in the license agreement ». Then click on « Next> ». A window appears :

Do not modify the access path on the above window, and click on the « Next> » button. On the following window click on « Install ».


If the below window appears during the installation, click on « OK » to continue.

Once the installation is finished, the window «Bluetooth favorites»All Bluetooth Neighborhood» appears. Simply close it. Then you see the following configuration window appear :

If you wish you can unselect the 3 boxes, and click on « Next> ». The below window appears :


Click on « Next> »

When the below window appears click on « Next> » again :

2 RSV-HD VIEWIRELESS

In the following window, uncheck all the service boxes excepted the service « **Bluetooth serial port** », then **click on « Next > »** :


In the below window **click on « Skip »** :



Finally, in the last installation window, **click on « Finish »** :

Now, we will setup the bluetooth peripheral. For this, click with the right button of the mouse on the Bluetooth icon in the notification zone. A menu is appearing, select the line « **Bluetooth configuration** ».



In the following window, **click on the « Client applications » tab** and select « **Bluetooth serial port** » from the list :

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
Click on the « **Properties** » button. The following screen appears :


Uncheck the box « **Secure connection** » and click on « **OK** ». Click on « **OK** » as well in the « Bluetooth configuration » window in order to close it.


A window named « My Bluetooth Places/Entire Bluetooth Neighborhood » appears. If not double click on the Bluetooth icon in the zone of notification.


Once you see the « **My Bluetooth Places** » window, start your RSV-HD Viewireless in picking up the sensor from its support (RSV-HD logo appears on the screen of the RSV-HD system). In the « **My Bluetooth Places** » window, click on the menu « **View** » and then « **Refresh** ». The computer detects the RSV-HD Viewireless system and it appears in the following way :

Then, double click on the icon with your unit's serial number. For example, « **RSV_VW_04056** ». The following window appears :

Double click on the icon « **Bluetooth the Serial Port on RSV_VW_04056** » in order to connect te computer with the RSV-HD Viewireless. The below window appears :

Windows mentions which serial port is used to connect the Bluetooth key. Please, remember you the serial port mentionned, you will need to use it in order to setup the peripheral in the Visident Imaging software. You can check (if you wish) the box « **Do not display this message again** », and click on « **OK** ».


2 RSV-HD VIEWIRELESS

Setup of RSV-HD Viewireless in Visident Imaging.


Start Visident Imaging in double clicking on the shortcut of your desktop :

In the imaging module, select the family « **Radiographs** » in the family list. You will use this family to work with your RSV-HD viewireless system.


In order to link your RSV-HD Viewireless system with Visident Imaging software, click on « **Parameters** », then on « **Acquisition** ». The following window appears :


In the type of acquisition list, select « **sensors** » by clicking on the « + » sign, then select the « **RSV-HD WIRELESS** » peripheral.

In this window, click as well on the « **Setup** » button. The following window appears :

Choice of the connection

The Wireless technology of the RSV-HD Viewireless system allows you to share your digital radiology system between several chairs. Select the configuration below which corresponds to your own use.

The sensor is used :

By a unique computer : Select this option if you want to use your RSV-HD Viewireless with a single computer. The other computers of the dental office will not be able to use the RSV-HD System, but it will be possible to share the images.

By computers in LAN : Select this option if you want to use the RSV-HD Viewireless with several computers. For this, computers have to be in a LAN. Thus, when a computer of the dental office needs to use the RSV-HD Viewireless, the user only needs to click on the activation button.

If another computer was already connected with the module, it is automatically disconnected to allow the new user connection.

By several computers not in LAN : Select this option if you want to drive the RSV-HD Viewireless with several computers. Because your computers are not in a LAN, you will have to manage manually the connection and the disconnection of the RSV-HD Viewireless module between each computer.

RSV-HD VIEWWIRELESS


Options

Connection timeout : This option is available only if you use the RSV-HD Viewwireless with several computers not in LAN (See choice of the connection above). If you choose to activate this option, the computer will automatically disconnect the RSV-HD system when the timeout over. Thus, another computer can connect the RSV-HD module without having to disconnect manually.

RSV-HD Bluetooth modules

You will define in the window below the RSV-HD Viewwireless module that you will use with your computer :


In order to add a module in the list, pick up the sensor from its support to start the activation of the system. Then, **click on « New »**.

The following window appears. Select the COM port in the list located on the top left of the window. This COM port has been provided to you at the end of the bluetooth key installation (see Chapter above). If you do not remember each of the COM Port, try one by one, that you have used, the COM Ports of the list.


Once the COM Port is selected, click on the button « Find the RSV-HD module ». Be sure that the RSV-HD Viewwireless module is still connected. When the peripheral is connected, you will see it appear in the list :


Then, click on the peripheral and **click on « OK »**.

2 RSV-HD VIEWIRELESS

A new window allows you to choose a label for this peripheral. Enter the label (or keep the one written by default) and **click on « OK »**.


The new peripheral is then appearing in the list of the wireless modules :

If you have other RSV-HD Viewireless modules in your dental office, please add them to the list using the same way than the first module.

Once all the modules are registered, select the module that you want to use by default with Visident Imaging. For this check the box of the peripheral you want :


Click on « OK » to close the configuration window.


Capture pictures

Now you are ready to capture pictures with your RSV-HD Viewireless.

Please note that the module you have chosen by default is selected in the family list :


The other modules that you have registered appear on this list as well. Just click on the one you want to use :


If you have setup your peripheral to be used with a unique computer, you are ready to take a picture. For this, pick up the sensor from its support. When the peripheral is synchronised with the computer, the orange LED appears on the RSV-HD Viewireless. You can now take a picture.

If you have setup your peripheral to use it with several computers, you must manually activate the connection with the RSV-HD Viewireless. For this click on the red activation button on the right of the family list.

When the activation button turns green your computer is connected with the RSV-HD module :

You are now ready to use RSV-HD Viewireless to take images.

For additional details on the use of Visident Imaging software, please refer to the manual provided on the CD.


- This RSV-HD user and installation manual

2 | 2 Use

2 | 2 | 1 Wireless connection with the computer

Step 1 : Sensor hung up


The system is in stand by mode, power supplier is off, only the sensor is supplied to detect the pick up of the CCD sensor.

② RSV-HD VIEWIRELESS


Step 2 : Sensor pick up

When you pick up the sensor of its location, the system activates thanks to an infra red cell and in two seconds the RSV-HD logo appears on the TFT screen for a duration of ten seconds. The blue LED blinks indicating the bluetooth connection with the PC (around 25 seconds). The orange LED is ON during a very short time as well (1 second) indicating that the internal card is correctly supplied.


Step 4 : Image processing

Once the capture has been made, image appears on the integrated TFT screen (after 3 - 4 seconds) for a duration of 30 seconds. For 6 to 10 seconds, the blue LED blinks quickly : it corresponds to the time to send the image in the PC via Bluetooth.

The « **one touch** » button is available only during this step. With short pressures, you can rotate the image and allows you to a long pressure (2 seconds) delete the image on the integrated screen and on the PC as well. During this long pressure, orange LED is on.

Step 5 : Back to step 3 or hang up the sensor and back to step 1

N.B. :

Orange LED is not blinking = the CCD sensor is not correctly plugged

Orange LED is blinking = Batteries are weak and need to be recharged.

2 | 2 | 2 Stand Alone Mode


The stand alone mode allows you to use the RSV-HD Viewwireless without any computer.

The pictures that you will obtain thanks to this mode will not be stored on the computer. However, you can do your diagnosis directly from the screen of the RSV-HD Viewwireless. A pressure on the « one touch » button will allow you to take a new picture.

The use is quite similar with the use of « wireless connection with the computer » described in the previous paragraph. Here is a description of the use in stand alone mode :

Step 1 : Sensor hung up


The system is in stand by mode, power supplier is off.

Step 2 : Activation in stand alone mode


To activate the system in stand alone mode :

- Press the « **one touch** » button and maintain the pressure.
- Pick up the sensor from its support (keep the pressure on the « one touch » button).

Step 3 : Initialization of the system


- During the initialization period, The orange LED is ON during a very short time as well (1 second) indicating that the internal card is correctly supplied.
- Then, logo appears on the screen.

② RSV-HD VIEWWIRELESS


Step 4 : End of the system initialization

Once the system is initialized (you can see the logo on the screen), stop the pressure on the « **One touch** » button.


Step 5 : Waiting for acquisition


Then, you can hear beeps indicating that the system is ready for the acquisition. In addition, the blue LED blinks slowly until the X-rays detection.

Step 6 : Display of the image

Once the capture has been made, image appears on the integrated TFT screen (after 3 - 4 seconds). Then it remains on the screen as long as you need.

Step 7 : Use of the « One Touch » button

You can use the « one touch » button in the following way :

- With short pressures, you can rotate the image
- With a long pressure, you delete the image

When the image has been deleted, you have two options :

If the sensor is not hung up, the system is initialized for a new capture (step 5).

If the sensor is hung up, the system returns in stand by mode (back to step 1).

N.B. :

Once you have finished using the system in stand alone mode, hang up the sensor and press for few seconds on the « one touch » button. The system returns stand by mode.

To use again, the system with the mode « wireless connection with the computer », you only have to pick up the sensor from its support (as described in the previous paragraph).

2 | 3 Batteries

Prior to the first use of your system RSV-HD Viewwireless or RSV-HD Wireless, it is strongly advised to charge the provided batteries (AAA batteries) as much as possible (until the charger LED is turned off).

Simply place the batteries in the charger while paying attention to the polarity. It is advised to use the providing charger (MASCOT 2116085000-131024-018110(EU)/2(UK)/4(Australia)/1(US)).

Moreover, after the first use of your system, be sure that the batteries are completely empty before recharging them.

Finally, it is advised to remove the batteries of the device if you are not using it for a long time.

The use of these batteries is without any risk and no particular care is needed.

214 Caution in the use of your sensor

The normal wear of your sensor can lead to damages. Thanks to the « hot plug & play » system (your sensor is automatically recognized by your system) you will not have to change to whole device in case of sensor wear.

RSV-HD sensors are high technology devices that need to be handled very carefully. Sensors have a 24 month warranty against eventual manufacturing or materials defects. Claims during the period of the warranty will be taken into account if the sensors have been used in normal conditions, and handled with care as mentioned below.

Caution :

The main cause of damage to CCD sensors is due to violent shocks, excessive pressure on the sensors, excessive bending or electrostatic discharges.

We remind you that the damages caused by this type of abuse are not covered by the warranty.

Caution: how to handle the sensor

In order to prevent damages that may cause problems with your sensor, read the following points :

What you can do :

- Wipe your sensor with a disinfecting towelette.
- Put back the sensor in a safe place when you are not using it
- Use a new hygienic protection cover with each use.

Do not :

- Do not sterilize the sensor with an autoclave
- Do not soak the sensor in a sterile solution
- Do not exert too much pressure on the sensor or its cable

- Do not pull the sensor by its cable when you want to remove the hygienic protection cover or the holder
- Do not hang the sensor by its cable
- Do not let your patient bite the sensor

Static electricity

Static electricity may be a cause of major damage on medical devices. It can be the case when you are using the sensor in an environment with carpets and dry air.

Thus, in order to avoid any problem we advise you to :

- Be sure that your electrical installation is correctly wired
- Install an antistatic floor-covering
- Use antistatic cleaning products
- Use an air humidifier

Caution in the handling of the cable

In order to prevent damages that may cause problems to your cable, read the following points :

What you can do :

- Untwist the cable carefully in letting the connector hang while you keep the sensor in your hands
- Be sure that the cable is not twisted before using
- Hold the connector and not the cable when you want to disconnect the sensor

Do not :

- Do not let the cable lie on the floor
- Do not let the cable hang near a drawer
- Do not clamp the cable
- Do not twist the cable while you are using it
- Do not pull the cable to disconnect the sensor

Caution :

Never use latex fingerstalls as hygienic protection for the sensor, but only hygienic protection covers recommended by VISIODENT. Fingerstalls contain talcum powder and are too tight on the sensor. The use of these fingerstalls may lead to tearing of the connection between the cable and the sensor that will not be covered by the warranty.

2 | 5 Caution in the use of your controller box

The controller box requires no specific maintenance apart from the usual cleaning and disinfection. It can be cleaned with a soft cloth. When cleaning, be careful not to let any liquid seep into the instrument, as this could cause permanent damage to the controller and permanently harm the way it operates. It is necessary to store the controller in a dry place.

2 RSV-HD VIEWIRELESS

2.16 Technical characteristics

	RSV-HD VIEWIRELESS
DIMENSIONS	178 mm x 60 mm x 33 mm
BACKLIT TFT SCREEN	●
ONE TOUCH SYSTEM	●
VERY HIGH DEFINITION MAX RESOLUTION SIZE 1 SIZE 2	1500 X 1000 PIXELS 1700 X 1200 PIXELS
WIRELESS BLUETOOTH TECHNOLOGY	●
BLUETOOTH USB KEY FOR PC	●
USB PORT 2.0 / 1.1 + USB CABLE 5m	
SUPPORT	●
FIXATION STRAP	●
IMAGING SOFTWARE	●
CHARGER 3.5 MM PLUG	●
4 AAA BATTERIES 1.2V EACH, CAPACITY 1600 mA/H	●

For all models: CE conformity / medical and electromagnetic security IEC 601-1


Voltage : 4,8 V / Current : 1,2 A

2.17 Fixation

FIXATION

You can use the wall-mount support or attach it to your equipment using the strap supplied with the RSV-HD. It can also be installed on any flat surface.

By fixing it to the chair or the arm of the x-ray unit, your system is always ready to use. This gives you greater comfort with your patient.

2 | 8 Information concerning the RSV-HD Viewwireless

Maximal external dimensions :

Length 17.8cm x width 6cm x depth 3.2cm
(+1cm depth corresponding to the support of the sensor).

LCD screen dimensions :

Length 4.5cm x width 3.5cm

Transmission :

Bluetooth Connection 2.0 EDR

Provided with a Bluetooth key to connect to the computer

Display time of the image on the TFT screen :

3 to 4 seconds

Transmission time :

6 to 7 seconds (this is not the display time)

Distance :

Maximum 10 meters with a space without obstruction

Power supply :

4 AAA batteries (1.2V and 800mAh each)

Autonomy :

Approximately 50 images

CE conformity :

Directive 93/42/EEC

Standards :

NF EN IEC 60601-1 for electromedical devices.

NF EN IEC 60601-1-2 for electromedical devices.

3 RSV-HD MEMORY


3.1 Installation and configuration

The installation procedure for this model is identical to that of the RSV-HD Viewwireless model in paragraph 2/1. Please refer to this chapter for the installation of your RSV-HD Memory model.


RSV-HD Memory Configuration in Visident Imaging

In the imaging module, select the family « **Radiographs** » in the family list. You will use this family to work with your RSV-HD memory system.


In order to link your RSV-HD memory system with Visident Imaging software, click on « **Parameters** », then on « **Acquisition** ». The following window appears :


In the type of acquisition list, select « **sensors** » by clicking on the « + » sign, then select the « **RSV-HD WIRELESS** » peripheral.

In this window, click as well on the « **Setup** » button. The following window appears :

Check the box « **RSV-HD Viewwireless Memory** ».

The wireless technology of RSV-HD Memory offers you the possibility to share your radiology system between different chairs in the dental office. Select the connection which corresponds to your use.


The sensor is used :

By one single computer : Select this option if you wish to use your RSV-HD Memory system with only one computer in the dental office. The other computers in the office will not be able to drive the RSV-HD system. However, they will be capable of sharing the patient's images taken with it.

By several computers in a network : Select this option if you wish to run your RSV-HD Memory system by several computers in the office. For this option, it is imperative that all of the computers in the dental office be linked in a network. When a computer in the dental practice is chosen to use the RSV-HD peripheral, the user must click on the activation button to make the connection. If a different computer has already established connection with the module, it will be disconnected automatically so that the new connection can be made.

By several computers not in a network : Select this option if you wish to operate your RSV-HD Memory system using several computers within the dental office. Since your computers are not linked in a network, you must manage the activation and deactivation of your RSV-HD Memory manually between the different computers.


Options

Connection timeout : This option is only available if you use the peripheral RSV-HD « with several computers not connected by a network » (see this option above). If you choose to activate this option, the computer automatically frees the Bluetooth peripheral once the selection delay is met. This way, another computer in the dental office can connect with the peripheral without first having to deactivate from another location.

Wireless peripherals

It is here that you will define the RSV-HD Memory peripheral(s) that you will use with your computer.


To add a module to the list, pick up the sensor from its support so that it will activate. Then, click on the button « New ».


The following window will appear. Select the Port COM from the list situated at the top left of the window. This Port COM has been furnished to you at the end of the installation of the Bluetooth key (see previous chapter). If you cannot remember which port you have selected, try each port on the list, one at a time.

Once the port has been selected, click on the button « Refresh the list » whilst assuring that the RSV-HD Memory peripheral it still activated. (To activate, hang up the sensor on its support and then pick it up again). When the peripheral is detected, it will appear on the list.


Select the peripheral by clicking on it and then clicking OK.

A new window allows you to choose a label for this peripheral. Enter the label (or keep the one written by default) and click on « OK ».


The new peripheral is then appearing in the list of the wireless modules :

If your dental practice has other RSV-HD Memory peripherals that you wish to use with your computer, add them to the list in the same manner as indicated for the first peripheral. Once all of the peripherals have been registered, select the peripheral that you would like to use by default in Visident Imaging. To do this, simply check the desired peripheral:


3 RSV-HD MEMORY


Click on « OK » to close the configuration window.

Taking images

You are now ready to take x-rays with your RSV-HD Memory peripheral. Notice that the RSV-HD Memory peripheral that you have selected as default is selected in the family list of images :


The other RSV-HD Memory peripherals should also appear when you scroll down the list. Simply click on the one which you wish to use.

If you have configured your peripheral for use with one single computer, you are ready to take an image. For this, pick up the sensor from its support. When the peripheral is synchronised with the computer, the orange light will appear on the controller box. You can now take the x-ray.

If you have configured your peripheral for use with several computers, you must manually activate the connection with the RSV-HD Memory peripheral. To do this, simply click on the red button situated next to the families list :

When this button turns green, your computer is connected with the RSV-HD Memory peripheral.


You are now ready to use the RSV-HD Memory peripheral to take your images. For more details on the use of the Visident Imaging software, Please refer to the Visident Imaging user's manual provided on the CD.

3 | 2 Use

3 | 2 | 1 Wireless connection mode with the computer


Step 1 : Sensor hung up

The system is in stand by mode, power supplier is off, only the sensor is supplied to detect the pick up of the CCD sensor.


Step 2 : Sensor Pick-up

When you pick up the sensor from its support, the system activates itself within two seconds thanks to the internal light detector. The RSV-HD logo should then appear on the monitor for a duration of ten seconds. The message « connection in progress » appears as well as the N° of the firmware of your module. Lastly, there is a battery indicator that informs you of the apparatus' charge. The blue LED blinks to signify the initialisation of the connection between the Bluetooth and the PC. (around 25 seconds). The orange LED will also light up for one second to indicate the card has a supply of energy.


**Step 3 : Wait for capture**

When the connection has been established, the Bluetooth symbol will appear on the module's monitor with the message « PC Mode ».

Then, you can hear beeps indicating that the system is ready for the acquisition. Moreover, the blue LED blinks slowly until the X-rays detection.

**Step 4 : Image treating**

Once the x-ray has been taken, the image will appear on the integrated monitor after 3-4 seconds. It will remain until either the sensor has been hung up or until the next image has been taken. The blue LED will flash rapidly for 6 to 10 seconds during the time that the image is being sent to the PC via Bluetooth.

3 RSV-HD MEMORY

The « **one touch** » button is available only during this step. With short pressures, you can rotate the image and will allow you to a long pressure (2 seconds) delete the image on the integrated screen and on the PC as well. During this long pressure, orange LED is on.

Step 5 : Back to step 3 or hang up the sensor and back to step 1

N.B. :

Orange LED is not blinking = the CCD sensor is not correctly plugged

Orange LED is blinking = Batteries are weak and need to be reloaded


3|2|2 Use in stand alone mode

The stand alone mode permits you to use the RSV-HD Memory module without a computer. The images taken in this mode are stocked in the internal memory of the module and can be eventually transferred to the computer when you wish to do so. The RSV-HD Memory in stand alone mode can be used in two ways :

3|2|2|1 Classic Stand Alone Mode (Without a PC nearby)

Step 1 : Sensor hung up


The system is in stand by mode, power supplier is off.

Step 2 : Activation in stand alone mode


To activate the system in autonomous mode :

- Press the « **one touch** » button and maintain the pressure.
- Pick up the sensor from its support (keep the pressure on the « one touch » button).

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- Whilst in the initialisation phase, the orange LED will quickly light up for 1 second to signify that the card has a supply of energy.
- The logo will then appear on the integrated monitor along with an indication of the level of battery power.

**Step 4 : End of controller box initialisation**

Once the controller box has been initialised, (logo and battery charge indicator showing on the screen) release the pressure on the button. The screen will show « Memory ».


Release the pressure on the button

This initialisation screen is followed by a new screen indicating the name of the dossier (represented by a letter of the alphabet) wherein the images of the current session will be stored.

Step 5 : Waiting for acquisition

Then, you can hear beeps indicating that the system is ready for the acquisition. Moreover, the blue LED blinks slowly until the X-rays detection

3 RSV-HD MEMORY

Step 6 : Image display

Two seconds after the x-ray has been taken, the image will appear on the integrated monitor. It will remain on the display for as long as you wish.

Return to step 5 if you wish to take additional x-rays. While this step is in process, you have the ability to see the x-rays taken by pressing (briefly) on the button in the centre of the module two times. The message « Display » will then appear on the module's screen.

In the Display mode, you can see all of the x-rays in the session by quickly pressing on the same button.

Pressing the button for a longer period will allow you to return to step 5. Hanging up the sensor on its support will end your current session.

Note :

When you restart your RSV-HD Memory module in classic stand alone mode, the last dossier will appear by default. If you wish to create a new dossier, simply press the button briefly when the dossier appears.

Attention :

As soon as the beep resounds, it is no longer possible to create a new dossier. In this case, please restart the manipulation.


3|2|2|2 Advance Stand Alone Mode (With a PC nearby)

In advanced mode, the files once labeled with letters of the alphabet (A, B, C, ...) are labeled with the names of your patients instead.

To send the patient's name to the module via the PC, it is necessary that the Bluetooth connection is established and that your Visident Imaging software is opened to the desired patient.

Step 1 : Sensor hung up

The system is in stand by mode, power supplier is off, only the sensor is supplied to detect the pick up of the CCD sensor.

Step 2 : Sensor pick-up

When you pick up the sensor from its support, the system activates itself automatically after 1 to 2 seconds thanks to a light detector. The RSV-HD logo will appear on the screen for approximately 10 seconds. The message « connection in progress » will appear with the N° of the firmware for your module. Finally, a battery indicator will inform you of the level of the system's charge. The blue LED will blink to signify that the Bluetooth connection with the PC is being initialised. (around 25 seconds) The orange LED will quickly light up for one second to indicate that the card has a supply of energy.

Phase 3 : Wait for capture


Once the connection has been established, the Bluetooth symbol will appear along with the message « PC Mode ».

Then, you can hear beeps indicating that the system is ready for the acquisition. Moreover, the blue LED blinks slowly until the X-rays detection.

TFT screen turned off


Blue LED blinks slowly until the x-rays detection

You can hear beeps

Sensor picked up

As soon as the beep sounds, you can see the following button appear on the Visident Imaging screen.


By clicking on the button, the name of the patient which you have opened is sent to the RSV-HD MEMORY system. The system moves into MEMORY mode (with the name of the patient displayed in red) and the connection with your computer is suspended.

The button  appears. You can now take as many x-rays as you would like for this patient.


If, for any reason, (for example if you wish to send a new patient name to the RSV-HD MEMORY system, if you would prefer to work in stand alone Memory mode or if you wish to work in PC mode) you wish to re-establish the Bluetooth connection with your computer, click again on the following button : . Refer to the different steps of the Classic Stand Alone mode.

3 | 2 | 2 | 3 Transfer mode

When you wish, you have the possibility of transferring your x-rays to your computer by clicking on the following button :

The following window will appear :

3 | 3 Batteries

Cf. rubrique 2 | 3

3 | 4 Precautions for sensor use

Cf. rubrique 2 | 4

3 | 5 Precautions for the use of the controller box

Cf. rubrique 2 | 5

3 | 6 Technical characteristics

TECHNICAL SPECIFICATIONS

DIMENSIONS	178 mm x 60 mm x 33 mm
BACKLIT TFT SCREEN	●
ONE TOUCH SYSTEM	●
VERY HIGH DEFINITION MAX RESOLUTION	Size 1 : 1500 x 1000 px Size 2 : 1700 x 1200 px
WIRELESS BLUETOOTH TECHNOLOGY	●
BLUETOOTH USB KEY FOR PC	●
SUPPORT HOLDER	●
FIXATION STRAP	●
IMAGING SOFTWARE	●
CHARGER 3.5 MM PLUG	●
4 AAA BATTERIES 1.2V EACH, CAPACITY 1,600 mA/H	●

For all models: CE conformity / medical and electromagnetic security IEC 601-1

3 | 7 Fixation

Cf. rubrique 2 | 7

3 | 8 Statistics of the RSV-HD Memory

Memory: Internal Micro card SD 1 Go

Maximal external dimensions :

Length 17.8cm x width 6cm x depth 3.2cm
(+1cm depth corresponding to the support of

LCD screen dimensions :

Length 4.5cm x width 3.5cm

Transmission :

Bluetooth Connection 2.0 EDR

Provided with a Bluetooth key to connect to the computer

Display time of the image on the LCD screen :

3 to 4 seconds

Transmission time :

6 to 7 seconds (this is not the display time)

Distance :

Maximum 10 meters with a space without obstruction

Power supply :

4 AAA batteries (1.2V and 800mAh each)

Autonomy :

Approximately 50 images

CE conformity :

Directive 93/42/EEC

Standards :

NF EN IEC 60601-1 for electromedical devices.

NF EN IEC 60601-1-2 for electromedical device.

4 RSV-HD WIRELESS


4 | 1 Installation and configuration

The installation and configuration process of this model is similar to the one describes for the RSV-HD Viewireless in paragraph 2.1. Please refer to this chapter for the installation and the configuration of your RSV-HD Wireless.

4 | 2 Use


Step 1 : Sensor hung up

The system is in stand by mode, power supplier is off, only the sensor is supplied to detect the pick up of the CCD sensor.


Step 2 : Sensor pick up

When you pick up the sensor of his location, the system activates thanks to an infra red cell. The blue LED blinks indicating the bluetooth connection with the PC (around 25 seconds). The orange LED is on for a short period of time (1 second) indicating that the internal card is correctly supplied.

Step 3 : Waiting for acquisition

Once the PC is connected (around 25 seconds) and the arming mode is activated (driven by Visident Imaging : the imaging module will appear automatically if you were previously on your patient file). Then, you can hear beeps indicating that the system is ready for the acquisition. Moreover, the blue LED blinks slowly until the X-rays detection. Once the PC is connected (around 25 seconds) and the arming mode is activated (driven by Visident Imaging : the imaging module will appear automatically if you were previously on your patient file), the blue LED blinks slowly until the X-rays detection.

Step 4 : Image processing


Once the capture has been made, for 6 to 10 seconds the blue LED blinks quickly : it is the necessary time to send the image to the PC via Bluetooth.

Step 5 : Back to step 3 or hang up the sensor and back to step 1

N.B. :

Orange LED is not blinking = the CCD sensor is not correctly plugged

Orange LED is blinking = Batteries are weak and need to be reloaded

**4 | 3 Batteries**

See information provided on [2 | 3](#)

**4 | 4 Caution in the use of your sensor**

See information provided on [2 | 4](#)

**4 | 5 Caution in the use of your controller box**

See information provided on [2 | 5](#)

4 | 6 Technical characteristics

	RSV ^{HD} WIRELESS
DIMENSIONS	178 mm x 60 mm x 33 mm
BACKLIT TFT SCREEN	
ONE TOUCH SYSTEM	
VERY HIGH DEFINITION MAX RESOLUTION	1500 X 1000 PIXELS
TABLE 1	1700 X 1200 PIXELS
TABLE 2	
WIRELESS BLUETOOTH TECHNOLOGY	●
BLUETOOTH USB KEY FOR PC	●
USB PORT 2.0 / 1.1 + USB CABLE 5m	
SUPPORT	●
FIXATION STRAP	●
IMAGING SOFTWARE	●
CHARGER 3.5 MM PLUG	●
4 AAA BATTERIES 1,2V EACH, CAPACITY 1 600 mA/H	●

For all models: CE conformity / medical and electromagnetic security IEC 601-1

Voltage : 4,8 V / Current : 1,2 A

4 RSV-HD WIRELESS

4 | 7 Fixation

See information in 2 | 7

4 | 8 Information concerning the RSV-HD Wireless

Maximal external dimensions :

Length 17.8cm x width 6cm x depth 3.2cm
(+ 1cm depth corresponding the support of the sensor).

Transmission :

Bluetooth Connection 2.0 EDR

Provided with a Bluetooth key to connect to the computer

Transmission time :

6 to 7 seconds (this is not the display time)

Distance :

Maximum 10 meters with a space without obstruction

Power supply :

4 AAA batteries (1.2V and 800mAh each)

Autonomy :

Approximately 50 images

CE conformity :

Directive 93/42/EEC

Standards :

NF EN IEC 60601-1 for electromedical devices.

NF EN IEC 60601-1-2 for electromedical devices.

5 RSV-HD USB

5 | 1 Installation

Your RSV-HD System is provided with its own «Visident Imaging» software». However, it is also possible to use this system with other imaging softwares thanks to some bridging tools that can be provided by Visident upon demand. Visident takes no responsibility for the operation of software other than Visident Imaging.

How to install Visident Imaging :


In order to use the RSV-HD peripheral you have to install the latest version of the Visident Imaging software. Please use the installation CD and follow the instructions in the installation manual of the software.

How to install the RSV-HD USB :


Please launch Visident Imaging by double clicking on the shortcut located on the desktop or on external imaging in Visident :


In the imaging module, select the « Radios » category in the family list. You will use your RSV-HD USB system with this family.

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
Please select « RSV-HD USB » and click on the Configuration button; the following window appears.

Click on « Install Driver ». Once it is installed, please click on « yes » to restart computer to complete install.

CAUTION

For the optimal use of your RSV-HD USB system, you must be sure that the electrical input of your computer is not disturbed with tension variations or micro cut. These problems will lead to a USB connexion loss. If you experience this problem, we advise the use of an uninterruptible power supply.


When the computer starts please connect your RSV-HD USB sensor.

5 | 2 Configuration


In order to configure your RSV-HD USB in Visident Imaging, please run Visident Imaging by double clicking on the shortcut located on the desktop or on external imaging in Visident :


Make sure the « Radios » family is selected.

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

When the following window appears please select « RSV-HD USB ».

Click on Configuration and this time click on « Auto detect COM ».

When the sensor is detected click on OK to close the window and click once again on OK to valid. Now you have the following window with a message « System ready » on the top of the screen.

5 RSV-HD USB

5 | 3 Use

Your sensor is installed and ready to use.

Step 1: Sensor hung up :

The system is in standby mode; the green LED of the controller box is off.


Location of the sensor

Step 2 : Sensor picked up :

When you are ready to take an exposure, start the system by clicking on the green switch on button on your software. The system triggers after 1 or 2 seconds and the green LED blinks of the computer ; thus the message « system Ready » appears on the screen.

Pick up of the sensor

Sensor picked up

Then you take radiographies.

The system automatically reloads after each acquisition, and is in standby mode when you hang up the sensor.

5 | 4 Caution in the use of your sensor

See information provided on 2/4

5 | 5 Caution in the use of your controller box

See information provided on 2/5

5 | 6 Technical characteristics

DIMENSIONS	150 mm x 50 mm x 22 mm
BACKLIT TFT SCREEN	
ONE TOUCH SYSTEM	
HD+ VERY HIGH DEFINITION MAX RESOLUTION SIZE 1	1500x1000 PIXELS
SIZE 2	1700x1200 PIXELS
WIRELESS BLUETOOTH TECHNOLOGY	
BLUETOOTH USB KEY FOR PC	
USB PORT 2.0 / 1.1 + USB CABLE 5m	
SUPPORT	
FIXATION STRAP	
IMAGING SOFTWARE	
CHARGER 3.5 MM PLUG	
4 AAA BATTERIES 1.2V EACH, CAPACITY 1 600 mA/H	

For all models: CE conformity / medical and electromagnetic security IEC 601-1

6 TECHNICAL CHARACTERISTICS OF THE CCD SENSOR

SENSORS S1/S2

ACTUAL RESOLUTION	> 20 pl/mm
THEORETICAL RESOLUTION	26,3 pl/mm
GREY SCALE	4096
PIXEL SIZE	19 μm^2
TECHNOLOGY	CCD + fiberoptics + CSI scintillator
NORMS	CE conformity / medical and electromagnetic security IEC 601-1

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USE	SIZE 1 periapical / others	SIZE 2 retroalveolar / bitewing
ACTIVE SURFACE AREA	30 mm x 20 mm (600 mm ²)	34 mm x 24 mm (816 mm ²)
DIMENSIONS	26.3 mm x 38.7 mm x 4.9 mm	42.6 mm x 30.6 mm x 4.9 mm
LENGHT OF THE SENSOR WIRE	80 cm	80 cm

The high sensitivity of these sensors makes it possible to reduce the dose of radiation. Supplied with a complete set of Visident sensor holders.

7 KIT OF POSITIONNERS

The digital radiology system of the RSV-HD line that you have bought is provided with a kit of positionners. These accessories are necessary to maintain the sensor at a right angle with the X-ray tube. They come in a separate case with their own user's manual. There are several kits fitted to the size of the sensor that you have. The positionners sold with your system must be sterilized with an autoclave.

8 WARRANTY

The systems of the RSV-HD line are proposed with a 24-months warranty at the date of invoice

This warranty does not cover :

- Any defects or damage caused by:
 - Shipment
 - Misuse of the product (i.e. the product is (i) used in inappropriate temperature or humidity or is exposed to excessive dirt, dust or fumes, (ii) sterilised with UV – light or (iii) cleaned with inappropriate solvents).
 - Water infiltrations into the system
 - Use of instruments incompatible with the product
 - Sudden mains power cuts or circumstances beyond our control such as fire or water damage or normal wear and tear...
- Any product that has been subject to abuse or negligence, to which any modifications, alteration or attachments have been made without our prior written approval or which has been installed or operated violating instructions for safe use.
- Any defect arising out of materials provided by, or a design stipulated or specified by the buyer.
- Any other expenses whatsoever, any personal injury or damage to any property, loss of production or profit or any other consequential or incidental economic loss.
- Any other equipment not manufactured by us.
Such equipment may carry a warranty given by the manufacturer thereof and if such warranty has been given to us, we will assign the benefit thereof to the buyer without recourse to us.

However, if the product has some material or manufacturing defects, Visident will repair this system and if it is not possible, Visident will proceed to its exchange. The customer will have to take in charge the shipping costs and the shipping insurance. Visident will take in charge the shipping costs and the shipping insurance back to the buyer.

NB : In the frame of the application of our warranty, the complete system must be shipped back in its original packaging.

Prescription only

This device is a prescription device and will only be used in healthcare facilities (e.g., dental offices) by trained personnel.

TO CONTACT US

Visiodent

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92210 La Plaine Saint-Denis FRANCE

Phone / + 33 (0) 1 49 46 58 00

Fax / + 33 (0) 1 49 46 58 08

Website / www.visiodent.com

Email / contactus@visiodent.com

FRANÇAIS

ENGLISH

Email / contact@visiodent.com

www.visiodent.com

VISIODENT