# IMAGING SOFTWARE

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Email / [contact@visiodent.com](mailto:contact@visiodent.com)

[www.visiodent.com](http://www.visiodent.com)

VISIODENT

# IMAGING SOFTWARE

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+

FRANÇAIS

## SOMMAIRE

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## WELCOME

You have recently purchased a digital radiology system from Visident group.  
We thank you for your trust and we will do our best to make you satisfy.

This user manual describes the way to use your Visident imaging software.

It contains all the information you will need on the software starting, the setup and the use of the imaging module.

Visident company shall not be held responsible for the consequences resulting from non respect of all the precautions describes in the present manual.

Visident, Visident XP, RSV - Radiology System Visident -, RSV HD, Visident Imaging, VPX, Elegance, Visiocam Wireless, as well as logos are Trademarks of the Visident Group.

Microsoft, Windows are Trademarks of the Microsoft Corporation.

All the other product names mentioned in the present manual are used for identification purposes and may be brands or trademarks registered by their owners.

This installation manual is the intellectual property of the Visident Group with all rights reserved. Any partial or complete reproduction of this manual in any form is strictly prohibited without the prior written consent of the Visident Group. Because of the constant evolution of its products, Visident reserves the right to modify the contents of this manual even without prior notice.

## 1 GETTING STARTED

Install Visiodent Imaging software with the installation manual provided with the CD. **Start the Visiodent Imaging in clicking on this icon.**

Once the Praticte and Practitioner parameters are entered, this window appears to access the Patient file.

Patient

New


### 1 | 1 Patient file

#### Creation of the patient file

- **Patient** : Seeking a Patient file
- **New** : Creation of a Patient file

#### Creation mode

Click on New to create a new patient file :

Imagerie


Enter : the Surname, and First Name, this is the minimum obligatory information to create a patient file (the rest can be entered later). **Click on the Floppy disk** at the top left to save.

To access the imaging module, **click on the imaging icon** on the right-hand bar.

## 2 IMAGING MODULE


### The menu

Located under the title bar, the menu proposes different headings.

### 2 | 1 File

« File » is composed of 4 sub menus.

### Configuration of print


Enables you to configure the printing to know how the plate(s) will be printed. (For more details, refer to paragraph III./ F/ a./ix.).

### Print view

Is used to visualise how the plate will come out on paper board before printing it.

## Print


Is used to print the plate presently displayed on the Visident Imaging body. A click on this obtains the following print window.

Select the printer on which you want to print and give the number of copies and **click on OK** to launch printing.

## 2 | 2 Parameters

« **Parameters** » is composed of the following sub menus.


### 2 | 2 | 1 Acquisition

This sub menu allows to parameter sensors for acquisition. By a simple click on the left button of the mouse, the following window appears :

## 2 IMAGING MODULE

### 2 | 2 | 1 | 1 RSV HD WIRELESS

#### Configuration of the RSV Viewwireless or Wireless module

The setup of the module is the same than the other peripherals. Select RSV HD Wireless in the “sensor in mouth” heading (see the screen above).


## 2 | 2 | 1 | 2 RSV Sensor

You can choose also the format to store the image :

In order to setup your RSV HD Viewless or Wireless peripheral, please refer to the manual provided with the peripheral.

You can choose also the format to store the image :

## 2 IMAGING MODULE

Then, click on « Configuration ». The following window appears :

The remaining time (default value : 10 mn) corresponds to the time the RSV System is connected.

If you check the box "hide the exposure information", the information indicating the time exposure will not appear on the image.

The box "High resolution acquisition" has to be selected if you want the default X-rays to be taken in high resolution. If not, by default it will be in low resolution.

### Details of the acquisition resolutions

#### High resolution


|                           |                |
|---------------------------|----------------|
| Image definition          | Very good      |
| Requested exposition time | 120 to 300 ms  |
| Acquisition time          | ~8 seconds     |
| Size of the image files   | ~100Ko in Jpeg |
| Quality with x2 zoom      | Very good      |
| Quality with x4 zoom      | Very good      |

#### Low resolution

|                            |               |
|----------------------------|---------------|
| Image definition           | Good          |
| Requested acquisition time | 40 to 120 ms  |
| Acquisition time           | ~2 seconds    |
| Size of the image files    | ~30Ko in JPEG |
| Quality with x2 zoom       | Very good     |
| Quality with x4 zoom       | Good          |

### Filters :

It enables to choose the filter to apply automatically during acquisition. Depending on the quality of operation of the Timer, Generator pair, the Gamma value will be documented by the installer when he has carried out the quality calibration tests.

**Recommended configuration :**  
Filter - Level Adjustment  
Gamma Level - 100

**Test**


Enables you to test the chosen configuration.

**Information**

Gives the information described below on the version of the equipment, the driver and various information about the product.

**2 | 2 | 1 | 3 Video for Windows**

According to the type of source chosen, the setup window is not the same. If you are in Video for Windows and that you **click on "setup"**, the following window appears :


This window is used to configure the video acquisition source. Here in the scroll list under the title bar, we choose the video capture driver to use. Once the driver is selected the desired source and format have to be defined.

Note : Screens copy of the following examples may be different according to the acquisition card installed in your computer.

Example :

**Click on Source**

Here you must give the type fo video connector (Composite or S-Video) and set the display properties (brightness, contrast, saturation), as well as the video standard used. (default value is: PAL B,D,G,H).

## ② IMAGING MODULE

It now remains to give the display standard. Once the source is established, you have to give the format in which the captures will be saved (dimension, colours).

Value to declare :

For image dimensions : Full

For image format : 24 bit RGB


Example :


Select **“overlay”** for image fluidity. The box “adjust image to window” is used to adjust the image to the body of Visident Imaging.

## 2 | 2 | 1 | 4 Direct show

« Direct Show » is currently the standard to drive imaging peripherals like video camera. This standard replaces progressively “Video For Windows”. To drive a camera, **choose Direct Show** and **click on setup**, the following window appears :


By clicking on **Crossbar video**, you choose the video entry to use. Check the type of connection of your camera, and choose the corresponding entry. Ex : « Video composite In » or « S-Video ».

By clicking on OK, you come back to the main setup window. Then, **click on Norme Video** to choose the kind of entry signal via the "standard video" list :

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
By clicking on OK, you come back to the main setup window. Then, **click on capture format** to choose the display format. Visiodent recommends to do the setup as you can see on the screen copy below:



## 2 IMAGING MODULE

### 2 | 2 | 1 | 5 Twain

Twain standard allows to drive peripherals as scanners, and many other peripherals. In order to use a TWAIN peripheral in Visident Imaging, click on TWAIN in the list, then **click on "setup"**. The following window appears :

**Caution :** Parameters available may be different according to the acquisition card installed on your computer.

Once you have chosen the peripheral, click on "Select", then click on "OK" to finish the setup.



## 2 | 2 | 1 | 6 VPX300 Panoramic & VPX400 Pano / Cephalo

These two Visident panoramics do not need setup. Select only the device in the list.

## 2 | 2 | 2 Display options

The display options are very useful to configure the display quality, using the buttons for screen setting, contrast and brightness, especially on flat screens and the screens of laptop computers (TFT, XGA, and so on...), this test is composed of 6 different widows which will enable you to determine in 3 tests (reference) if the angle of your screen is correct for optimum use of your RSV Imaging. Tests are :


Reference 1


If on the left image the same shapes can be distinguished as on the right (background, angles and edges), then the angle is correct.

## 2 IMAGING MODULE



Reference 2

This test consists in doing the same as the previous test. Gamut

Reference 3


In the right part of the window, you can distinguish very clearly the black square in the middle of the yellow as well as the white square in the middle of the blue. If you can distinguish these same 2 squares in the window on the left with the same sloping angle, then it is correct. Different levels of grey can also be distinguished, from the darkest to the lightest all around the window.

## Status bar

## 2 | 2 | 3 Colours selection

By clicking on this sub menu, the following window appears :

This window allows to customized colours used in this software :  
Background colour of the images. Background colour of the windows.

## 2 | 2 | 4 Convert the images

If you select this sub menu, the images stored in the format TIFF and JPEG from previous versions of Visident Imaging will be converted in new PNG and JPEG2000 format.

## 2 | 2 | 5 Status bar


This option allows to display or to hide the status bar of Visident Imaging. This bar provides useful information during the use of Imaging as the number of available images, the brightness, contrast, gamma values and so on...

## 2 | 2 | 6 Palettes

Select here the palettes to use in order to display x-rays in pseudo colours. The choice of the palettes can be done thanks to a right click on the pseudo colours icon as well.

## 2 | 2 | 7 Zoom configuration

Allows to modify the zoom characteristics :

You can modify :




- The zoom factor from 2 to 6
- The size of the magnifying glass from 32 to 128

The configuration window of the zoom is available as well by a right click on the zoom icon.

## 2 IMAGING MODULE

### 2 | 2 | 8 Flashlight configuration


Allows to modify the flashlight characteristics. You can determine the size and the intensity of the filter.

| Types  | Setup  | Détails  |
|--|--|--|
| <input checked="" type="radio"/> Brightness/contrast mode<br><input type="radio"/> VISIODENT(r) filter<br><input type="radio"/> Gamma mode<br><input type="radio"/> Histogram equalization | <input checked="" type="radio"/> Brightness/contrast mode<br><input type="radio"/> VISIODENT(r) filter<br><input type="radio"/> Gamma mode<br><input type="radio"/> Histogram equalization<br>Flashlight size<br> 96<br>Brightness value from -128 to 128<br> -32<br>Contrast value from -128 to 128<br> 32 | Allows to determine the size of the flashlight in mm. The 2 other tools give the intensity, the brightness and the contrast. |
| <input checked="" type="radio"/> VISIODENT(r) filter   | Flashlight size<br> 96  | To know the size of the flashlight in mm.  |
| <input checked="" type="radio"/> Gamma mode  | Flashlight size<br> 96<br>Gamma value (x100) from 0 to 500<br> 200   | To know the size of the flashlight in mm and the gamma value.  |
| <input checked="" type="radio"/> Histogram equalization  | Flashlight size<br> 96  | Allow to know the size of the flashlight in mm.  |

The flashlight setup window is available by a right click on the mouse on the flashlight icon too.

### 2 | 2 | 9 Rendering configuration

Allow to choose how the image must be displayed on the screen :

**Adjusted display** : displays the images as large as possible in full display.

**1:1 Display** : displays the image with its real size.


**Optimised display** : displays the image at the size of the best rendering.

The rendering configuration window is available with a right click on the mouse above the display images zone too. (on the middle of the screen).


## 2 | 2 | 10 Status configuration

This menu is available only in "dental status display" mode. This window allows to choose in a list the model of dental status with which you want to work. You have to choose as well the options you want to display in your status.

## 2 | 2 | 11 3D display configuration

This menu is only available in 3-D display mode. This window allows to select the 3-D display that you wish :



## 2 IMAGING MODULE

### 2 | 2 | 12 Organisation of the favourite tab setup


This sub-menu allows you to choose until 12 functions that will appear on the favourite tab. These functions can be management, tools or processing functions.

In order to add a function in this tab:

1. Select the function to add (a red frame appears on the tool)
2. Click on the pointer in order to move the icon from the left to the right
3. Then icon appears on the favourite tab


Use the pointers located on the bottom of the screen to move the function where you want :

You can add and delete easily the functions with a « **drag & drop** ». For this :


1. Select an icon in maintaining the button of the mouse pushed.
2. Move the icon to the place you want. The place appears in a red frame.
3. Release the button to fix the icon on the chosen place.

When you have finished to setup the favourite tab, you just have to click on the cross to close the window :


Then, functions of the favourite tab can be used as you can do for the other tabs management, tools or processing.

## 2 IMAGING MODULE


### 2 | 3 Language

Here choose the language that you want to use in Visident Imaging :

### 2 | 4 Menu ?

« ? » menu give you information on Visident Imaging.


#### Help

This sub-menu enables to obtain necessary information about the way to use your software. Moreover, the "F1" key enables to open the contextual help. When you position the mouse on a specific heading and you press "F1" key, corresponding information of this heading appear on the screen.

#### About Visident Imaging

By clicking on this sub menu, you obtain the following window. This window gives information on Visident Imaging: (version, ...).

If you have an Internet access, you can send an e-mail to Visident team in order to submit all your suggestions on future software evolutions.

## 2 | 5 Acquisition bar

Description of the acquisition bar

### Family list

(or combo box) located in the upper right side of the window, allows to know in which image family you are.

If you want to work on radiographs or to capture images, you have to select "**Radiographs**".


If you want to work on colour images or do a film to capture images, select "**Photographs**".

For Panoramic radiographs, please select "**Panoramics**".

If you want to work on all the images, radiographs, colour images or panoramics, select <ALL> family. Be careful, you cannot do any capture with this family. <TOUTES> en sachant que cette famille ne permet aucune acquisition.

If your practice has a RSV HD Wireless system (model Viewwireless or Wireless), each recorded peripheral will appear in the family list.

In the below list, the practice has 2 RSV HD Wireless peripheral that you can see in the family list for each computer of the practice.


Note :To setup your RSV HD Wireless, please refer to the manual provided with the system.


The user can then switch easily from one to the other peripheral. For this, you only have to select the peripheral with which the computer has to connect.

## ② IMAGING MODULE

If the user wishes to use the peripheral  
**RSVHD WIRELESS 1 :**


Now if the user needs to use the peripheral  
**RSVHD WIRELESS 2 :**


### Acquisition buttons for the RSV HD Wireless peripheral


The use of your peripheral depends of the type of connexion that you use. You can know it in going in menu **"Parameters"** -> **"Acquisition"** and then in clicking on **"setup"**:


« **By a unique computer** » : For this connexion mode, there is no acquisition button, because the connexion between the computer and the RSV HD Wireless is completely automatic.

« **By computer in LAN** » : In this mode, you need to activate manually the connexion with the RSV HD Wireless.. For this you need to click on the red button located close to the family list :

When the button changes for green light, it means that your computer is connected with the RSV HD Wireless peripheral :

## ② IMAGING MODULE

- ① If the peripheral is already connected to another computer, this connexion is automatically cut and the green button becomes red (it means that the connexion with the peripheral is over). The computer « A » uses the RSV HD Wireless peripheral :


- ② When the user of the computer « B » wants to use the RSV HD Wireless peripheral, he clicks on the activation button. It becomes green. Computer « A » is automatically disconnected. The activation button becomes red :

- ① **"By computer not in LAN":** With this mode, computers cannot communicate together, users must activate and deactivate manually the connexion with the RSV HD Wireless. Let's take an example of the use of one RSV HD Wireless system with two computers in this mode. Computer "A" use the RSV HD Wireless:

- ② Computer « A » do not need anymore to use the peripheral and disconnects it in pressing the activation button. Button becomes red. RSV HD Wireless is not connected to any computer :

- ③ Once ce peripheral disconnected, computer « B » can use the peripheral. For this, the user need to click on the activation button. Button becomes green.


For additional details on the RSV HD Viewwireless or Wireless, please refer to the manual provided with the device.

## 2 IMAGING MODULE

### Acquisition buttons for the RSV peripheral

In radio family, to drive the RSV sensor, available buttons are :

The first button allows to connect or to disconnect the sensor. The two following buttons allow to switch between high and low resolution.

### Acquisition buttons for WDM or DirectShow

In Photographs family, you can drive a camera with a Direct Show (WDM) format. Here are the available buttons :

With the first button you can play or stop the video stream coming from the camera. With the second button you can take a picture. The captured image is then stored in the patient file.

### Acquisition buttons for a panoramic

In Panoramics family, to drive a Panoramic VERAVIEWPOCS, VPX300 or VPX400 :

This button allows to activate the panoramic or cephalometric device according to the model you have.

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Fonctions bar of Visident Imaging is composed of 4 buttons: Favourite, Management, Tools and Processing. Click on one of the four buttons to see the corresponding functions displayed.

### Favourite buttons

After the software installation, this button contains no functions. You have the possibility to setup 12 functions in this tab. Functions available are the one you can find in the other tabs : Management, Tools and Processing. Once your functions have been setup in the Favourite tab, you can find them easily without to look for them in several tabs. To setup the functions, please refer the Favourite part above.

Management

Import

Export

Image  
properties

### Management buttons

#### Import

Is used to import an image into Visident Imaging. This button is the only one present even if no image is selected.

#### Export

Is used to export images out of the imaging directory of Visident Imaging.

#### Image properties

To make comments on an image, to give it tooth/teeth number(s) and to give a date to the image.

Portrait

**Portrait**

Take a picture of your patient and click on this button to see the picture displayed in the patient file of this patient.

Quotation

**Quotation**

Only available with Visiodent Imaging version linked with Visiodent management software. By clicking on this button, Visiodent Imaging do a quotation on the selected image if this quotation is setup in the management software.

Cancel

**Cancel**

Is used to cancel the process to be carried out (one by one). This button is available in the three buttons of the fonctions bar.

Cancel all

**Cancel all**

Is used to cancel all the processes carried out since the acquisition of the image.

Delete

**Delete**

Is used to delete an image.

Print

**Print**

Is used to print an image. You can modify the printing configuration in the menu "File", then "print configuration"; or by clicking with the right button of the mouse on the icon.


Preview

**Preview**

Is used to preview the printing rendering before to print it definitively.

By clicking on the right button of the mouse on the preview icon, you obtain the following screen:

With this window you can configure the printing to know the display of the pictures before to print.

## 2 IMAGING MODULE

Tools

### Tools buttons

Rotation -90°

#### Rotation - 90°

Is used to turn the image 90° to the left.

Rotation +90°

#### Rotation + 90°

Is used to turn the image 90° to the right.

Rotation +90°

#### Rotation 180°

Is used to turn the image through 180°.

Horizontal flip

#### Miroir horizontal

Is used to turn the image around a horizontal axis.

Vertical flip

#### Miroir vertical

Is used to turn the image around a vertical axis.


Histogram

### Histogram

The histogram represents the number of pixel depending on the light intensity, from the black on the left, to the white on the right. It is used to know if the image had a good exposure. The graph should have a curve like a bell and the top of this graph should be in the middle of the image has a good exposure.

The tools “fit to width of window” and “centre of window” are used to modify the histogram, that is to say the brightness in the image. Thus we can give contrast to images without contrast and reduce the contrast on image with too many contrast.

Zoom

**Zoom**

Is used to magnify a chosen part of the image :


The zoom parameters can be set. To do this, make a right click on the zoom icon. You can use the sub menu **"zoom configuration"** in the **"Parameters"** menu as well.

Flashlight

**Flashlight**

Is used to apply a chosen filter on an image.


The flashlight size as well as the filters can be set. A right click with the mouse on the flashlight allows to open the configuration window. You can use the sub menu **"Flashlight configuration"** in **"Parameters"** menu as well.

## 2 IMAGING MODULE

### Isocontour extraction


Is used to highlight image zones corresponding to a specific grey level.

### Length measurement

By clicking on this button, you will be able to carry out a measurement on the image. You must remember that "the measurements are for your information only and should not be used for diagnostic purposes."

By clicking on the image at the first time, you have a departure point for the measurement, then by clicking a second time you designate the end. The measurements appear in the right frame and the sum at the bottom.

### Angle

Is used to measure the angle between two segments. Click on the extremity of both segments to see the angle between these two segments.

### Processing buttons

Section

### Section

With this button by drawing a line on the image, you will be able to determine which elements are at the same grey level.

Negative

### Negative

Is used to reverse the values of the image.

Bright. Contrast  
Gamma

### Brightness contrast Gamma

By clicking on this icon, you make appear brightness, contrast and Gamma tools in the cursor bar along the fonctions bar. The same bar is used for pseudo colours and for isocontours extraction functions.

LN Correction

### LN Correction

Is used to read a too white image.

Smoothing

### Smoothing

Is used to smooth in case of under exposure.

Pseudo  
colours

### Pseudo colours

Is used to highlight the grey levels in colour. This filter is useful when grey levels are very similar. Replacing these grey levels enables to better distinguish the undertones of an images.

Periodontics  
filter

### Periodontics filter

Is used to highlight the gum on the images.

Endodontics  
filter

### Endodontics filter

Is used to highlight the roots on the images.

Dentine filter

### Dentine

Is used to highlight the crown on the images.

Visiodent filter

### Visiodent filter

Is used to highlight several relief levels.

Niveaux auto

### Level adjustment


Give contrast to the image to highlight details.

## 2 IMAGING MODULE


### Normal display : one image

This button is used to switch the software over into normal display (on single image). Example :

### Comparison mode (two images)

This button is used to switch the software over into comparison display (two images). Example:

In this mode the left frame moves to the top and you simply click on an image in this frame to display it.


- A right click and the image appears in the frame on the right
- A left click and the image appears in the frame on the left.

### Dental status display

This button is used to switch the software over into dental status mode. Example :

New status

The dental status is composed of five specific buttons specific :

**New status** : start a new status, all the frames are empty.

Status archiving

**Status archiving** : Is used to record images as well as their positions. The status displays the date of the backup and appears on the list below the action buttons.

Delete status

**Delete status** : Is used to delete the selected status in the status list.

Automatic filling

**Automatic filling** : Is used to fill in the status automatically. For this you must previously have saved the images with teeth numbers.

Personalised status

**Personalised status** : Under construction.

Empty the base

**Empty the base** : Is used to delete an image in a status. However, this image will not be delete in the patient's imaging file.


**Backup of the status** : All the modifications that you bring to a status are automatically backup as soon as you leave the status mode. The saved status mentioned the date allowing to find it easily later.

## ② IMAGING MODULE


### Full screen display

This button is used to switch the software over into full screen display mode. Example :

### Patient's image file display


Example :

### 3D display

This button is used to switch the display over into 3D mode.

### Exit

The left frame displays the saved images of a current patient according to the selected family.

### The left frame

In this frame, images are displayed as well as 3 buttons :

#### Buttons

|                                     |                 |
|-------------------------------------|-----------------|
| <input checked="" type="checkbox"/> | Sort by date    |
| <input type="checkbox"/>            | Sort by name    |
| <input type="checkbox"/>            | Ascending sort  |
| <input checked="" type="checkbox"/> | Descending sort |

To choose the sort criterion for display.

Under construction.

To make the left frame disappear or reappear to enlarge the display of the body frame.

## ② IMAGING MODULE

### 2 | 7 Images

Below these 3 buttons, the images appear :


If you make a left click on it, the image is framed in yellow and appears in the Visident Imaging body frame or in the left part in comparison mode.

If you make a right click, the image is framed in blue and appears in the right part only in comparison mode.

### 2 | 8 The body

The Visident Imaging body is used to display a selected image, for this you have to click with the left click on the mouse and it appears in the Visident Imaging body.

## TO CONTACT US

Visiodent

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92210 La Plaine Saint-Denis

Phone / + 33 (0) 1 49 46 58 00

Fax / + 33 (0) 1 49 46 58 08

Website / [www.visiodent.com](http://www.visiodent.com)

Email / [contactus@visiodent.com](mailto:contactus@visiodent.com)

Email / [contact@visiodent.com](mailto:contact@visiodent.com)

[www.visiodent.com](http://www.visiodent.com)

VISIODENT